

Eye on HIV Testing

EVOLUTION & PRACTICAL ISSUES

UNAIDS / WHO HIV testing strategies*

Which Strategy to Follow Depends upon Prevalence*

Objective of testing	Prevalence of infection	Strategy applicable
Screening of blood & blood products	-	I
Surveillance	> 10	I
	< 10	II
Diagnosis	> 30	I
	< 30	II
Asymptomatic	> 10	II
	< 10	III

How to report HIV test results?*

Negative

If initial/screening test shows non-reactive result.

Positive

If the sample shows reactive results concordantly by the three screening tests.

Indeterminate (equivocal)


If the sample shows discordant results by the 3 screening tests, the follow up samples are required to retest at 2 weeks and at 3, 6 and 12 months before the final status of the test results should be conveyed. If the status remains indeterminate after 1 year, the person is considered to be HIV antibody negative.

Special Considerations

Pre-test, post-test counseling service should be provided and confidentiality should be maintained.


STRATEGY-I

Transfusion safety
Transplant safety
Surveillance




STRATEGY-II

Surveillance
Diagnosis



STRATEGY-III

Diagnosis



1. Assay A1, A2 and A3 represent three different assays.
2. Such a result is not adequate for diagnostic purposes; use strategy II or III. Whatever the diagnosis, donations which were earlier reactive should not be used for transfusion or transplant.
3. Report; result may be reported.
4. For newly diagnosed individuals a positive result should be confirmed on a second specimen.
5. Testing should be repeated on a second specimen taken after 14 days.
6. Result is considered negative in the absence of any risk of HIV infection.

Duration of window period

Authorised body	Duration of window period
CDC, USA	2 weeks-1 year; Average: 3-6 months
The Canadian Medical Association	Within 6 months, 95% seroconvert within 3 months
WHO	6 weeks-3 months, unusual cases: 6 months or longer
NACO, India	3 weeks-3 months on an average, Can be longer sometimes

Generations & Detectability of HIV Assays

HIV assay	Antigen used	Detects	Can Detect by
1st Generation	Viral lysate	IgG antibody	6-12 weeks after exposure
2nd Generation	Recombinant antigens	IgG antibody	6-12 weeks after exposure
3rd Generation	Recombinant antigens / Synthetic peptides	Total HIV antibody	3-4 weeks after exposure
4th Generation	Recombinant antigens / Synthetic peptides / mono &/or polyclonal p24 antibody	Total HIV antibody & p24 antigen	2-3 weeks after exposure
p24 antigen immunoassay	Mono &/or polyclonal p24 antibody	p24 antigen	2 weeks after exposure
Western Blot	Viral lysate/recombinant	IgG antibody	Eq. to 1st/2nd generation
Nucleic Acid Technology (NAT)	-	HIV RNA	1 week after exposure

Causes of false-positive result

Flu, Herpes simplex, Upper respiratory tract infection, Recent viral infection, Exposure to viral vaccines, pregnancy in multiparous women, Malaria, Rheumatoid arthritis, Hepatitis B vaccination, Tetanus vaccination, Autoimmune diseases, Haemolyzed serum, Tuberculosis etc.

Causes of false-negative result

- **Performance and technical errors :**
Pipetting error, mislabeling of samples, variability in test kits, powder from powdered gloves etc.
- **Biologic, pathologic, and pharmacologic determinants :**
Window (pre-seroconversion) period, delayed antibody synthesis in infants, immunosuppressive therapy, Concurrent infection with EBV /CMV etc.
- **Congenital or drug-induced hypogammaglobulinemia :**
- **Sensitivity and specificity of the assay:**
Sampling prior to immunoglobulin M to immunoglobulin G class switching (only if 2nd generation assay is used for testing), haemodilution, limited antigenic determinants etc.

Factors influencing the duration of window period

- The host-virus interaction
- Immune status of the host
- Route of infection
- Type / subtype of infecting HIV strain
- Viral inoculum at the time of exposure

Immunoassays for HIV Testing

QUALISA-HIV
ELISA to detect HIV 1 & 2 antibodies in human serum/plasma

RETROLISA 3.0
3rd generation ELISA to detect HIV 1 & 2 antibodies in human serum/plasma

RETROQUIC
Rapid immunoconcentration assay to detect & differentiate HIV 1 & 2 antibodies in human serum/plasma

RETROCHECK
3rd generation rapid ICT to detect HIV 1 & 2 antibodies in human serum/plasma/whole blood

RETROSCREEN
3rd generation ICT to detect & differentiate HIV 1 & 2 antibodies in human serum/plasma

PRACTICAL SOLUTIONS TO PRACTICAL ISSUES



QUALPRO DIAGNOSTICS



email: sales@tulipgroup.com
website: www.tulipgroup.com

*Strategy I is recommended for blood donation purposes.
 For the use of a registered medical practitioner or a hospital or a laboratory only. Refer package insert before use.