

PERFORMANCE CHARACTERISTICS

A) Internal Evaluation:

- Accuracy: In an internal study **Electra™ LH** was evaluated against commercially available licensed kit with 90 random clinical samples & **Electra™ LH** has demonstrated >98% clinical correlation with the commercially available licensed kit.
- Precision: **Electra™ LH** was evaluated with licensed external Quality controls for Precision Studies & following is the data:

Controls	No. of testings	Mean Control values with Electra™ LH	Coefficient of Variable (CV)
Level 1	10	1.488	8.61
Level 2	10	25.18	5.70
Level 3	10	92.58	4.44

B) External Evaluation:

Electra™ LH CLIA has been evaluated by a NABL accredited lab against their reference method. In this evaluation **Electra™ LH** has demonstrated 98% correlation with the reference method.

*Data file: Zephyr Biomedicals (A Division of Tulip Diagnostics Pvt. Ltd).

Important Note:

- The wash procedure is critical. Insufficient washing will result in poor precision and falsely elevated absorbance readings.
- It is recommended to use the multi channel pipettes to avoid time effect. A full plate of 96 wells may be used if automated pipette is available.
- Duplication of Standards & Samples is not mandatory but may provide information on reproducibility & application errors.












LIMITATIONS OF THE ASSAY

(1). As with all diagnostic tests, a definite clinical diagnosis should not be based on the results of a single test, but should only be made by the physician after all clinical and laboratory findings have been evaluated. (2). The activity of the enzyme used is temperature-dependent and the RLU values may vary. The higher the room temperature (+18°C to +25°C) during substrate incubation, the greater will be the RLU values. Corresponding variations apply also to the incubation times. However, the Standards are subject to the same influences, with the result that such variations will be largely compensated in the calculation of the result. (3). Adaptation of this assay for use with automated sample processors and other liquid handling devices, in whole or in part, may yield differences in test results from those obtained using the manual procedure. It is the responsibility of each laboratory to validate that their automated procedure yields test results within acceptable limits. (4). Insufficient washing (e.g., less than 5 wash cycles, too small wash buffer volumes, or shortened reaction times) can lead to incorrect results.

BIBLIOGRAPHY

- Knobil, E. The neuroendocrine control of the menstrual cycle, Rec. Prog.Horm. Res. 36: 52-88; 1980
- Harris, G.W. and Naftolin. The hypothalamus and control of ovulation. Brit. Med. Bullet. 26: 1-9; 1970
- Shome, B. and Parlow, A.F. J. Clin. Endocrinol. Metab. 39:199-202; 1974
- Shome, B. and Parlow, A.F. J. Clin. Endocrinol. Metab. 39:203-205; 1974
- Uotila, M.; Ruoslahti, E. and Engvall, E. J. Immunol. Methods. 42: 11-15; 1981

SYMBOL KEYS

 Temperature Limitation	 Consult Instructions for use	 Date of Manufacture	 Batch Number / Lot Number
 Manufacturer	 In vitro Diagnostic Medical Device	 This side up	 Contains sufficient for <n> tests
 Use by	 Catalogue Number	 Do not reuse	

Manufactured by:
Zephyr Biomedicals

A Division of Tulip Diagnostics (P) Ltd.

M 46-47, Phase III B, Verna Industrial Estate, Verna, Goa- 403 722, INDIA.

Regd. Office: Gitanjali, Tulip Block, Dr. Antonio Do Rego Bagh, Alto Santacruz, Bambolim Complex P.O., Goa - 403 202, INDIA.

An ISO 13485
Certified Company

0320/VER-01



Chemiluminescence Assay for the Quantitative Determination of Luteinizing Hormone (LH) in Human Serum.
FOR IN VITRO DIAGNOSTIC USE ONLY
Store at 2°C to 8°C

INTENDED USE

ELECTRA™ LH CLIA test is intended for the quantitative determination of Luteinizing Hormone (LH) in human serum. For In Vitro Diagnostic Use only.

INTRODUCTION

Luteinizing hormone (LH) is produced in both men and women from the anterior pituitary gland in response to luteinizing hormone-releasing hormone (LH-RH or Gn-RH), which is released by the hypothalamus. The alpha chain is similar to that found in human thyroid-stimulating hormone (TSH), follicle-stimulating hormone (FSH) and human chorionic gonadotropin (hCG). The differences between these hormones lie in the amino acid composition of their beta subunits, which account for their immunological differentiation.

The basal secretion of LH in men is episodic and has the primary function of stimulating the interstitial cells (Leydig cells) to produce testosterone. The variation in LH concentrations in women is subject to the complex ovulatory cycle of healthy menstruating women and depends on a sequence of hormonal events along the gonado-hypothalamic-pituitary axis. In the differential diagnosis of hypothalamic, pituitary, or gonadal dysfunction, assays of LH concentration are routinely performed in conjugation with FSH assays since their roles are closely interrelated. Furthermore, the hormone levels are used to determine menopause, pinpoint ovulation, and monitor endocrine therapy.

PRINCIPLE

ELECTRA™ LH Quantitative CLIA assay is for use on **ELECTRA** analyzers. **ELECTRA™ LH CLIA** works on the principle of chemiluminescence wherein light is produced by a chemical reaction from a substance as it returns from an electronically excited state to the ground state. When catalysed by HRP, the oxidation of luminol by hydrogen peroxide produces an electronically excited form of 3-aminophthalate which on relaxation emits light with maximum intensity at $\lambda=425\text{nm}$.

The **ELECTRA™ LH** Quantitative Test Kit is based on a solid phase enzyme immunoassay. The assay system utilizes one anti-LH antibody for solid phase (microtiter wells) immobilization and another mouse monoclonal anti-LH antibody in the antibody-enzyme (horseradish peroxidase) conjugate solution. The test sample is allowed to react simultaneously with the antibodies, resulting in the LH molecules being sandwiched between the solid phase and enzyme-linked antibodies. After incubation, the wells are washed with Wash Buffer to remove unbound labeled antibodies. A solution of chemiluminescent substrate is then added and Luminescence is measured in RLU. The intensity of the emitting light is proportional to the amount of enzyme present and is directly related to the amount of LH in the sample. By reference to a series of LH standards assayed in the same way, the concentration of LH in the unknown sample is quantified.

MATERIALS & COMPONENTS

Materials provided with the test kits:

- Coated Microwells: Microwells coated with Anti- LH antibody.
- LH HRPO Enzyme Conjugate
- Substrate A: Chemiluminescent substrate containing enhanced luminol solution.
- Substrate B: Chemiluminescent substrate containing stabilized peroxide solution.
- LH Standard set of 6 standards labeled as A to F in lyophilized form. Ready to use. For standard Concentrations refer vial label.
- Wash Buffer Concentrate (20X).

Materials required but not provided:

- Precision pipettes: 10-100 μl , 20-200 μl , 100-1000 μl
- Disposable pipette tips
- Distilled water
- Disposable Gloves
- ELECTRA™ Analyzer**

STORAGE AND STABILITY

- ELECTRA™ LH** kit is stable at 2-8°C up to the expiry date printed on the label.
- Coated micro-wells should be used within one month of opening the pouch. Once opened, the pouch must be sealed properly to protect from moisture. In case the desiccant pouch changes color from blue to white, the strips should not be used.
- Diluted wash buffer is stable up to one week at 2-8°C.
- Working Substrate (A+B) must be used immediately.



Chemiluminescence assay



Chemiluminescence assay



Chemiluminescence assay



Chemiluminescence assay

SPECIMEN COLLECTION

1. Collect blood specimen by venipuncture according to the standard procedure.
2. Only serum should be used.
3. Avoid grossly hemolytic, lipemic or turbid samples.
4. Preferably use fresh samples. However, specimens can be stored up to 48 hours at 2-8°C, for short duration.
5. For longer storage, specimens can be frozen at -20°C. Thawed samples must be mixed prior to testing.
6. Do not heat inactivate before use.
7. Specimen containing precipitate or particulate matter should be clarified by centrifugation prior to use.
8. Specimen should be free from particulate matter and microbial contamination.

PRECAUTIONS

1. Bring all reagents and specimen to room temperature before use.
2. Do not pipette any material by mouth.
3. Do not eat, drink or smoke in the area where testing is done.
4. Use protective clothing and wear gloves when handling samples.
5. Use absorbent sheet to cover the working area.
6. Immediately clean up any spills with sodium hypochlorite.
7. All specimens and standards should be considered potentially infectious and discarded appropriately.
8. Neutralize acid containing waste before adding hypochlorite.
9. Do not use kit after the expiry date.
10. Do not mix components of one kit with another.
11. Always use new tip for each specimen and reagent.
12. Do not allow liquid from one well to mix with other wells.
13. Do not let the strips dry in between the steps.


REAGENT PREPARATION

- All reagents should be brought to room temperature (18-25°C) and mixed by gently inverting or swirling prior to use. Do not induce foaming.
- Dilute Wash Buffer 20 times (for example add 5 ml concentrated buffer to 95 ml distilled or deionized water). Mix well before use.
- Prepare a Working Substrate by Mixing Substrate A and Substrate B in equal volume (1:1 ratio) before addition to the micro-wells.
- Since the reference standards are lyophilized, reconstitute each standard with 0.5ml distilled water. Allow the reconstituted material to stand for at least 20 minutes. Reconstituted standards should be sealed and stored at 2-8°C.

No. of Strips	1	2	3	4	5	6	7	8	9	10	11	12
Substrate-A µl	250	450	650	850	1050	1250	1450	1650	1850	2050	2250	2450
Substrate-B µl	250	450	650	850	1050	1250	1450	1650	1850	2050	2250	2450

TEST PROCEDURE

1. Secure the desired number of coated wells in the holder. Dispense **25 µl** of Standards & Serums into the appropriate wells.
2. Dispense **50 µl** of Enzyme Conjugate into each well. Incubate at room temperature (18-25°C) for **30 minutes**.
3. Remove the incubation mixture by emptying the plate content into a waste container. Rinse and empty the microtiter plate 5 times with Wash Buffer (1X). Strike the microtiter plate sharply onto the absorbent paper or paper towels to remove all residual water droplets.
4. Add **50 µl** of working Substrate (A+B) in all the micro-wells. Keep away from direct light while adding the substrate.
5. Cover the **ELECTRA™** microplate and incubate for **10 minutes** at room temperature (18-25°C) in dark.
6. Read the **ELECTRA™** micro-plate exactly at **10 minutes** in **ELECTRA™ Analyzer**. If **ELECTRA™** micro-plate is not read between 10-15 minutes the test results should be considered as invalid.



CALCULATIONS


Construct a standard curve by plotting the RLU obtained from each reference standards against its Concentrations in mIU/ml on the graph paper, with RLU values on the vertical or Y axis and concentrations on the horizontal or X axis. Use the RLU values for each specimen to determine the corresponding concentration of LH in mIU/ml from the standard curve. Any diluted specimens must be corrected by the appropriate dilution factor.

Example of Standard curve

Results of a typical standard run with RLU's shown in the Y axis against LH concentrations shown in the X axis.

Suggest: Use 4-Parameter Standard curve to calculate sample values.

LH Values (mIU/ml)	RLU's
A	1268
B	6010
C	117662
D	531039
E	1143856
F	2361492



This Standard curve is for the purpose of illustration only and should not be used to calculate unknowns. Each user should obtain their own Standard curve and data.

Expected Ranges of values

Each laboratory must establish its own normal ranges based on patient population. The results provided below are based on randomly selected clinical laboratory samples.

LH (mIU/ml)				
	Age	No. of patients	Mean	Range
Male	< 10	25	1.3	< 2.5
Male	15-60	56	4.8	1 - 15
Female	< 10	25	1.1	< 2.0
Female	20 - 35	60	15	1 - 90
Female	46 - 60	40	38	8 - 120

The minimal detectable concentration of human LH by this assay is estimated to be 2 mIU/ml.