

STORAGE AND STABILITY

1. **ELECTRA™ Tp. Ab. CLIA** kit is stable at 2-8°C up to the expiry date printed on the label.
2. Coated micro-wells should be used within one month of opening the pouch. Once opened, the pouch must be sealed properly to protect from moisture. In case the desiccant pouch changes color from blue to white, the strips should not be used.
3. Diluted conjugate must be used immediately.
4. Diluted wash buffer should preferably be used within the same day. However excess wash buffer may be stored at up 2-8°C for one week.
5. Working Substrate (A+B) must be used immediately.

MATERIAL REQUIRED BUT NOT PROVIDED

1. Manual or automatic pipettor.
2. Pipettor tips.
3. Incubator.
4. Micro-well washer
5. **ELECTRA™ SA** or **ELECTRA™ FA**.
6. Reagent grade water.
7. Disposable gloves
8. Timer/ Stop Watch.

SAMPLE COLLECTION

1. No prior preparation of the patient is required.
2. Collect blood specimen by venipuncture according to the standard procedure.
3. Serum or plasma can be used.
4. Specimen should be free of particulate matter and microbial contamination.
5. Use of fresh sample is preferred. However, specimen can be stored refrigerated for upto three days. For long storage, samples should be frozen at 2°C or below.
6. Specimen should be brought to room temperature prior to testing.
7. Do not heat inactivate before use.

PRECAUTIONS

1. Do not pipette any material by mouth.
2. Do not eat, drink or smoke in the area where testing is done.
3. Use protective clothing and wear gloves when handling samples.
4. Immediately clean up any spills with sodium hypochlorite.
5. Dispose off all the reagents and material used as if they contain infectious agent.
6. Neutralize acid containing waste before adding hypochlorite.
7. Do not use the kit after the expiry date.
8. Do not mix the components of one kit with those from another.
9. Always use a new tip for each specimen and reagent.
10. Do not allow liquid from one well to mix with other wells.
11. Do not let the strips dry in between the steps.
12. Do not expose working substrate to direct light.

REAGENT PREPARATION

1. Dilute wash buffer 20 times (for example add 5 ml concentrated buffer to 95 ml distilled or deionized water).
2. Dilute enzyme conjugate with Conjugate diluent according to the requirement as shown below. Prepare a fresh dilution for each assay.

No. of Strips	1	2	3	4	5	6	7	8	9	10	11	12
Enzyme Conjugate (µl)	40	80	120	160	200	240	280	320	360	400	440	480
Conjugate Diluent (µl)	960	1920	2880	3840	4800	5750	6720	7680	8640	9600	10560	11520

3. Prepare a Working Substrate by Mixing Substrate A and Substrate B in equal volume (1:1 ratio) before addition to the micro-wells.

No. of Strips	1	2	3	4	5	6	7	8	9	10	11	12
Substrate –A (µl)	250	450	650	850	1050	1250	1450	1650	1850	2050	2250	2450
Substrate –B (µl)	250	450	650	850	1050	1250	1450	1650	1850	2050	2250	2450

TEST PROCEDURE

1. Bring all the reagents and specimen to room temperature before use.

2. Take out required number of strips and immediately close the pouch.
3. Well A1 must be used as blank.
4. Add 100 µl Diluted conjugate in each well except well A1
5. Add 50 µl each of NC in wells B1 & C1 & PC in wells D1 & E1 respectively.
6. Add 50 µl of samples to the remaining wells.
7. Apply plate sealer and incubate for 45 minutes at 37°C.
8. Wash each well using a semi-automated micro-plate washer (Preferably Lisa Wash models) with diluted wash buffer for 6 wash cycles giving 30 secs soak time for each wash cycle and Blot dry.
9. Add 50 µl of working Substrate (A+B) in all the micro-wells. Keep away from direct light while adding the substrate.
10. Cover the Electra microplate and incubate for 10 mins at room temperature (18-25°C) in dark.
11. Read the Electra micro-plate exactly at 10 mins in **ELECTRA™ SA** or **ELECTRA™ FA**. If Electra micro-plate is not read between 10- 15 mins, the test results should be considered as invalid.

TEST VALIDATION CRITERIA

1. The individual RLU of negative controls should be less than 200000 for **ELECTRA™ SA** and **ELECTRA™ FA**.
2. The individual RLU of positive controls should be more than 7500000 for **ELECTRA™ SA** and **ELECTRA™ FA**.

CALCULATIONS

The cutoff value is calculated by the software of **ELECTRA™ SA** and **ELECTRA™ FA** using the mean of negative control RLU and a factor.

Factor for **ELECTRA™ SA** and **ELECTRA™ FA** is 400000

SAMPLE DATA

For **ELECTRA™ SA** and **ELECTRA™ FA**:

Considering the Cutoff for this batch is 400527

Well	RLU	ECI	Interpretation
Blank	80		
NC	541		
NC	513		
PC	14401045		
PC	14562704		
Sample 1	3780160	9.43	Reactive
Sample 2	9936	0.02	Non-reactive
Sample 3	1487376	3.71	Reactive
Sample 4	209727	0.52	Non-reactive

INTERPRETATION OF RESULTS

1. Results are interpreted in E.C.I. This is determined by dividing the RLU of the sample by the Cutoff value calculated for that specific run.
2. Samples with E.C.I values greater than or equal to 1.00 are considered Reactive and samples with E.C.I values less than 1.00 are considered Non-reactive.
3. Samples that are initially reactive in **ELECTRA™ Tp. Ab. CLIA** should be retested in duplicate. Repeated reactivity is highly predictive Syphilis infection.

PERFORMANCE CHARACTERISTICS

1. 449 samples out of which four *T. pallidum* positive & 445 *T. pallidum* negative samples were tested with **ELECTRA™ Tp. Ab. CLIA** and with compared with commercially available 3rd generation *T. pallidum* EIA.

Specimen Data	Total	ELECTRA™ Tp. Ab. CLIA	<i>Other T. pallidum</i> EIA
Total Specimens	449	449	449
Syphilis Reactive	4	4	4
Syphilis Non-Reactive	445	445	445

2. **ELECTRA™ Tp. Ab. CLIA** was also evaluated with 105 samples from RPR Syphilis Positive Research Panels (Panels obtained from SLR research, Carlsbad, CA).
3. **ELECTRA™ Tp. Ab. CLIA** was evaluated using In House syphilis Positive Panels (P.P/SY/E-09, & P.P/S/E-10) for sensitivity using serial dilution against other *T. pallidum* EIA.