

Size : 274 x 218 mm



**Enzyme Linked Immunosorbent Assay (ELISA) for detection of antibodies to HIV 1 & 2  
in human serum or plasma.**

**FOR IN VITRO DIAGNOSTIC USE ONLY**

**Store at 2°C to 8°C**

#### **INTENDED USE**

**QUALISA™ HIV 1/2** is intended to be used for the detection of antibodies to HIV 1 / 2 virus in human serum or plasma.

#### **SUMMARY AND EXPLANATION**

**QUALISA™ HIV 1/2** is an Enzyme Linked Immunosorbent Assay (ELISA), which employs highly purified synthetic peptides representing most conserved antigenic segments of envelope glycoproteins, gp 120 and gp 41 of HIV-1 and gp 36 of HIV-2.

#### **PRINCIPLE OF THE ASSAY**

Microwell strips are coated with synthetic peptides gp 120, gp41 and gp36 representing both HIV-1 and HIV-2. Diluted samples along with positive and negative controls are added in the coated wells and incubated. The wells are washed to remove unbound components and Agglutinating sera for Human IgG conjugated to horseradish peroxidase (HRPO) is added. After a short incubation the wells are washed again and bound enzyme is detected by adding substrate. The reaction is stopped after specified time with acid and absorbance is determined for each well at 450 nm with an ELISA reader. The cut-off value (COV) is calculated by the given formula and absorbances of all the wells are compared with the COV. Any sample having absorbance more than the COV is considered reactive.

#### **KIT COMPONENTS**

**QUALISA™ HIV 1/2** has following components:

1. Coated microwells : Microwells coated with synthetic peptides representing both HIV 1 and HIV 2. Ready to use.  
96 Wells: (3x8) x 4 pouches  
192 Wells: (3x8) x 8 pouches  
480 Wells: (3x8) x 20 pouches
2. Positive control: Inactivated and stabilized human serum reactive for HIV-1 and non-reactive serum for HCV & HBsAg with preservatives.
3. Negative control: Inactivated and stabilized human serum non reactive for HIV-1 and HIV-2, HBsAg and HCV.
4. Conjugate: Agglutinating sera for Human IgG-HRPO conjugate (50 X). To be diluted 50 times with conjugate diluent.
5. Conjugate diluent: Buffered solution containing stabilizing proteins and preservatives. Ready to use.
6. Sample diluent: Buffered solution containing stabilizing proteins and preservatives. Ready to use.
7. Substrate: Solution containing Tetramethyl benzidine (TMB) and hydrogen peroxide. Ready to use..
8. Wash buffer: Buffer contains surfactants (20 X). To be diluted 20 times with distilled or deionized water
9. Stop solution: Diluted sulphuric acid. Ready to use.
10. Microwell holder.
11. Instruction for use.
12. ELISA protocol sheet.
13. Plate sealer.

#### **STORAGE AND STABILITY**

1. **QUALISA™ HIV 1/2** kit is stable at 2-8°C upto the expiry date printed on the label.
2. Coated microwells should be used within two months of opening the pouch. Each pouch contains 3 strips of 8 wells. Any unused well(s) should be resealed with cellotape (and desiccant pouch) to prevent moisture absorption and stored at 2-8°C for future use.
3. In case the desiccant pouch changes colour from blue to white, the strips should not be used.
4. Diluted conjugate must be used immediately.
5. Diluted wash buffer is stable upto one week.

#### **MATERIAL REQUIRED BUT NOT PROVIDED**

1. Manual or automatic pipette.
2. Pipette tips.
3. Glass or polypropylene container for conjugate dilution.
4. Absorbent sheets.
5. ELISA washer.
6. ELISA reader.
7. Pipetting troughs or boats.
8. Disinfectant.
9. Reagent grade water.
10. Disposable gloves
11. Timer
12. Biohazard waste container
13. Serological pipettes

Size : 137 x 218 mm

#### SAMPLE COLLECTION

1. No prior preparation of the patient is required.
2. Collect blood specimen by venipuncture according to the standard procedure.
3. Serum or plasma can be used.
4. Specimen should be free of particulate matter and microbial contamination.
5. Preferably use fresh sample. However, specimen can be stored refrigerated for short duration.  
For long storage, freeze at -20°C or below. Do not freeze samples in frost-free freezer.
6. Specimen should not be frozen and thawed repeatedly.
7. Do not heat inactivate before use.
8. Specimen containing precipitate or particulate matter should be clarified by centrifugation prior to use.

#### PRECAUTIONS

(1). Bring all reagents and specimen to room temperature before use. (2). Do not pipette any material by mouth. (3). Do not eat, drink or smoke in the area where testing is done. (4). Use protective clothing and wear gloves when handling samples. (5). Use absorbent sheet to cover the working area. (6). Immediately clean up any spills with sodium hypochlorite. (7). Dispose off all the reagents and material used as if they contain infectious agent. (8). Neutralize acid containing waste before adding hypochlorite. (9). Do not use kit after the expiration date. (10). Do not mix components of one kit with another. (11). Always use new tip for each specimen and reagent. (12). Do not let the dispensing tip of ELISA washer touch liquid in the wells. (13). Do not allow liquid from one well to mix with other wells. (14). Do not let the strips dry in between steps. (15). Incubation times and temperatures other than those specified may give erroneous results. (16). Cross contamination of reagents or samples may give false results. Samples must remain confined to microwells during testing. (17). Inadequate removal of residual wash buffer can cause inconsistent color development. Microwell strips should be tapped vigorously and blotted on absorbent paper or towels to minimize residual wash buffer. (18). Reusable glassware must be washed and thoroughly rinsed free of detergents prior to use. All glassware must be clean and dry before use. (19). Adequate washing of wells is extremely important. Inadequately washed wells will exhibit high background values, and may show false positive values. For manual washing, aspirate the contents of the wells, then fill each well with wash buffer solution. Avoid cross-contamination of the wells, particularly in the first wash after aspiration. Drain all of the wash buffer from the wells by inverting, then shaking residual wash buffer from the wells with a sharp "snapping" motion of the wrist. Repeat these steps for at least six washes. The wells should then be tapped vigorously on a paper towel or other absorbent material to remove all traces of residual wash buffer. The use of an automated microwell washing system will assure consistent washing of the wells, and is recommended.

#### REAGENT PREPARATION

1. Dilute wash buffer 20 times (for example add 5 ml concentrated buffer to 95 ml distilled or deionized water).
2. Dilute conjugate 50 times (for example add 20 µl concentrated conjugate to 980 µl conjugate diluent).


No. of Strips	1	2	3	4	5	6	7	8	9	10	11	12
50 x Conjugate	20µl	40µl	60µl	80µl	100µl	120µl	140µl	160µl	180µl	200µl	220µl	240µl
Conjugate Diluent	980µl	1960µl	2940µl	3920µl	4900µl	5880µl	6860µl	7840µl	8820µl	9800µl	10780µl	11760µl

#### TEST PROCEDURE

1. Bring all the reagents and specimen to room temperature before use.
2. Take out required number of strips and immediately close the pouch.
3. Prepare ELISA protocol sheet indicating the location of controls and specimen.
4. Use controls in duplicate.
5. Add **200 µl** sample diluent to separate wells.
6. Add **10 µl** of control or specimen in the wells.
7. Gently shake the plate to mix thoroughly. Care should be taken to avoid spillage.
8. Apply plate sealer and incubate for 30 minutes at 22-28°C.
9. Wash each well by filling approximately **350 µl** diluted wash buffer and aspirating/flicking off six times. Blot dry.
10. Add **100 µl** diluted conjugate in each well and incubate for 30 minutes at 22-28°C.
11. Wash six times as in step 9. Blot dry.
12. Add **100 µl** substrate and incubate at 22-28°C away from light for 30 minutes.
13. Stop reaction by adding **100 µl** stop solution. The stop solution should be added in the same sequence as substrate addition.
14. Read the absorbance at 450 nm with 600-700 nm as reference within 30 minutes of stopping the reaction.

#### RUN CRITERIA

1. The individual absorbance value of negative controls should be less than 0.1.
2. The individual absorbance value of positive controls should be more than 1.0.



### CALCULATIONS

The cut-off value (COV) is calculated by adding 0.2 to average absorbance value of negative control.

#### EXAMPLE

$$COV = Av.NC + 0.2$$

Absorbance of Negative Control (NC)	
Reading-1	0.025
Reading-2	0.031
Average NC reading	0.028

Absorbance of Positive Control (PC)	
Reading-1	2.000
Reading-2	2.200
Average PC reading	2.100

$$\text{Cutoff value} = \text{Av. NC} + 0.2 = 0.028 + 0.2 = 0.228$$

### SAMPLE DATA

Well	Absorbance	Mean	Cutoff	Result
NC	0.025	0.028	0.228	
NC	0.031			
PC	2.000	2.100		
PC	2.200			
Sample 1	3.487			Reactive
Sample 2	0.654			Reactive
Sample 3	0.056			Non-Reactive

### INTERPRETATION OF RESULTS

1. Samples with absorbance value less than the cut-off value are considered non-reactive by **QUALISA™ HIV 1/2** are considered negative for HIV 1/2 antibodies.
2. Samples with absorbance value equal to or greater than cut-off value are considered reactive by **QUALISA™ HIV 1/2**. The original sample should be retested in duplicate. Initially reactive sample that do not react in either of duplicates are considered negative for antibodies to HIV 1 / 2. Initially reactive sample that reacts in either or both duplicates are considered repeatedly reactive.
3. If a sample is repeatedly reactive the probability of antibodies to HIV 1 / 2 are high, specially with patients at high risk or high absorbance values. Such samples should be retested with supplemental tests western blot or HIV RNA. Specimens that are repeatedly reactive in ELISA but not reactive in additional testing are considered indeterminate and a further sample after 3 to 6 months should be tested.
4. In case of samples with high OD, there are possibilities of black precipitate formation after the addition of stop solution. This will not interfere with the interpretation of results.

### NOTE

All ELISA readers have a linear measuring range (approx. 2.5<sup>0</sup>A). Beyond this range the OD values are non-linear. That is why some instruments can be pre-programmed to show "out of range" values if the OD exceeds 2.5 or 3.00. **QUALISA™ HIV 1/2** positive control may show "out of range" values in such instruments even after dilution. Please note that this indicate a valid run (provided the negative control value is < 0.1) and does not interfere with assay results.

### PERFORMANCE CHARACTERISTICS

One thousand one hundred and twelve specimens – out of which one hundred and fourteen positive and nine hundred and ninety eight negative specimens were tested with **QUALISA™ HIV 1/2** and compared with a commercially available ELISA kit utilizing similar principle. The results are given below:

SPECIMEN DATA	TOTAL	QUALISA™ HIV 1/2	Commercially available ELISA
Number of specimen tested	1112	1112	1112
Number of Positive specimens	114	114	114
Number of Negative specimens	998	997	997

Based on this evaluation:

Sensitivity of **QUALISA™ HIV 1/2** : 100%

Specificity of **QUALISA™ HIV 1/2** : 99.9%

### Evaluation with Seroconversion Panel

Results of Seroconversion Panel obtained from Boston Biomedica Inc. (PRB904) are given below:

Panel ID	Days since First bleed	Abbott HIV 1/2	Gen. Sys HIV 1/2	Org. Tek. HIV	Dupont WB HIV RNA	Roche PCR	QUALISA™ HIV 1/2
PRB 904-01	0	0.1	0.2	0.5	NEGATIVE	NEGATIVE	0.25
PRB 904-02	21	0.2	0.2	0.4	NEGATIVE	NEGATIVE	0.25
PRB 904-03	49	0.2	0.3	0.5	NEGATIVE	POSITIVE	0.35
PRB 904-04	92	12.0	3.5	5.1	POSITIVE	POSITIVE	5.0
PRB 904-05	99	12.2	4.7	4.9	POSITIVE	POSITIVE	12.5

**Note:** All numerical values are expressed as specimen absorbance to cut off ratios. Ratios more than or equal to 1.0 are considered positive. The data other than that for **QUALISA™ HIV 1/2** was supplied by BBI.

### Precision

Inter & Intra-assay precision studies indicate that the CV is 10 ± 3%.

### REMARKS

(1). Though **QUALISA™ HIV 1/2** is a reliable screening assay, it should not be used as a sole criterion for diagnosis of HIV infection. (2). Absence of antibodies to HIV does not indicate that an individual is absolutely free of HIV infection as the collection of sample and its timing vis-à-vis seroconversion will influence the test outcome. (3). Since various tests for HIV differ in their performance characteristics and antigenic composition, their reactivity patterns may differ. (4). Testing of pooled samples is not recommended. (5). As with all diagnostic tests, a definitive clinical diagnosis should not be based on the result of a single test, but should only be made by the physician after all clinical and laboratory findings have been evaluated. (6). **QUALISA™ HIV 1/2** should only be used as a screening test and its results should be confirmed by other supplemental method before taking clinical decisions.

### BIBLIOGRAPHY

(1) Centers for Disease Control, Update on Acquired Immune Deficiency Syndrome (AIDS) MMWR 1982; 31:507-508. (2) Popovic. M., et.al., Detection Isolation and continuous production of Cytopathic Retroviruses (HTLV-III) from patients with AIDS and pre-AIDS. Science 1983;224:497-500. (3) Harshida .S. et.al., Earlier detection of Human immunodeficiency Virus Type 1 p24 antigen and immunoglobulin G and M antibodies to p17 antigen in seroconversion serum panels by immune complex transfer enzyme immunoassays, Clinical and Diagnostic Laboratory Immunology, Nov 2000 p872-881, Vol 7, No.6. (4) Dorn et.al., Analysis of Genetic Variability within the Immunodominant Epitopes of Envelope gp41 from Human Immunodeficiency Virus Type 1 (HIV-1) Group M and Its Impact on HIV-1 Antibody Detection, Journal of Clinical Microbiology, Feb 2000, p.773-780. Vol.38. No.2. (5) Vanhems P. et.al., HIV seroconversion interval and demographic characteristics ; no evidence of selection bias: Sex Transm inf. 2002;77:446-448. (6) Cecilia C. L. Ngan et.al., Alternative Strategies for Confirmation of Human Immunodeficiency Virus Infection Require Judicious Use, Journal of Clinical Microbiology, Jan 2002, p 314-315, Vol 40, No.1. (7) Persaud D et.al., Latency in Human Immunodeficiency Virus Type 1 Infection: No Easy Answers, Journal of Virology, Feb 2003, p. 1659-1665. Vol. 77, No.3. (8) Data on file: Qualpro Diagnostics.



Manufactured by:

**Qualpro Diagnostics**

A Division of Tulip Diagnostics (P) Ltd.

88/89, Phase II C, Verna Industrial Estate, Verna, Goa - 403 722, INDIA.

**Regd. Office:** Gitanjali, Tulip Block, Dr. Antonio Do Rego Bagh, Alto Santacruz, Bambolim Complex P.O., Goa - 403 202, INDIA.

0617VER-02