

2. Precision: **Qualisa™ LH** was evaluated with licensed external Quality controls for Precision Studies & following is the data:

Controls	No. of testings	Mean Control values with Qualisa™ LH	Coefficient of Variable (CV)
Level 1	10	1.488	8.61
Level 2	10	25.18	5.70
Level 3	10	92.58	4.44

B) External Evaluation:

Qualisa™ LH ELISA has been evaluated by a NABL accredited lab against their reference method. In this evaluation **Qualisa™ LH** ELISA has demonstrated 98% correlation with the reference method.

*Data file: Zephyr Biomedicals (A Division of Tulip Diagnostics Pvt. Ltd.).

IMPORTANT NOTE

1. The wash procedure is critical. Insufficient washing will result in poor precision and falsely elevated absorbance readings.
2. It is recommended to use the multi channel pipettes to avoid time effect. A full plate of 96 wells may be used if automated pipetting is available.
3. Duplication of standards & samples is not mandatory but may provide information on reproducibility & application errors.












LIMITATIONS OF THE ASSAY

1. As with all diagnostic tests, a definite clinical diagnosis should not be based on the results of a single test, but should only be made by the physician after all clinical and laboratory findings have been evaluated.
2. The activity of the enzyme used is temperature-dependent and the OD values may vary. The higher the room temperature (+18°C to +25°C) during substrate incubation, the greater will be the OD values. Corresponding variations apply also to the incubation times. However, the standards are subject to the same influences, with the result that such variations will be largely compensated in the calculation of the result.
3. Adaptation of this assay for use with automated sample processors and other liquid handling devices, in whole or in part, may yield differences in test results from those obtained using the manual procedure. It is the responsibility of each laboratory to validate that their automated procedure yields test results within acceptable limits.
4. Insufficient washing (e.g., less than 5 wash cycles, too small wash buffer volumes, or shortened reaction times) can lead to incorrect OD values.

BIBLIOGRAPHY

1. Knobil, E. The neuroendocrine control of the menstrual cycle, **Rec. Prog. Horm. Res.** 36: 52-88; 1980
2. Harris, G.W. and Naftolin. The hypothalamus and control of ovulation. **Brit. Med. Bulletin.** 26: 1-9; 1970
3. Shome, B. and Parlow, A.F. **J. Clin. Endocrinol. Metab.** 39:199-202; 1974
4. Shome, B. and Parlow, A.F. **J. Clin. Endocrinol. Metab.** 39:203-205; 1974
5. Uotila, M.; Ruoslahti, E. and Engvall, E. **J. Immunol. Methods.** 42: 11-15; 1981

SYMBOL KEYS

	Temperature Limitation		Consult Instructions for use		Date of Manufacture		Batch Number / Lot Number
	Manufacturer		In vitro Diagnostic Medical Device		This side up		Contains sufficient for <n> tests
	Use by		Catalogue Number		Do not reuse		



Manufactured by:

Zephyr Biomedicals

A Division of Tulip Diagnostics (P) Ltd.

M 46-47, Phase III B, Verna Industrial Estate, Verna, Goa - 403 722, INDIA.

Regd. Office: Gitanjali, Tulip Block, Dr. Antonio Do Rego Bagh, Alto Santacruz, Bambolim Complex P.O., Goa - 403 202, INDIA.



Luteinizing Hormone

Enzyme Linked Immunosorbent assay for the Quantitative Determination of Luteinizing Hormone (LH) in Human Serum

FOR IN VITRO DIAGNOSTIC USE ONLY

Store at 2°C to 8°C

INTENDED USE

Qualisa™ LH Sandwich ELISA test is intended for the quantitative determination of luteinizing hormone (LH) in human serum. For In Vitro Diagnostic Use only.

INTRODUCTION

Luteinizing hormone (LH) is produced in both men and women from the anterior pituitary gland in response to luteinizing hormone-releasing hormone (LH-RH or Gn-RH), which is released by the hypothalamus. The alpha chain is similar to that found in human thyroid-stimulating hormone (TSH), follicle-stimulating hormone (FSH) and human chorionic gonadotropin (hCG). The differences between these hormones lie in the amino acid composition of their beta subunits, which account for their immunological differentiation.

The basal secretion of LH in men is episodic and has the primary function of stimulating the interstitial cells (Leydig cells) to produce testosterone. The variation in LH concentrations in women is subject to the complex ovulatory cycle of healthy menstruating women and depends on a sequence of hormonal events along the gonado-hypothalamic-pituitary axis. In the differential diagnosis of hypothalamic, pituitary, or gonadal dysfunction, assays of LH concentration are routinely performed in conjugation with FSH assays since their roles are closely interrelated. Furthermore, the hormone levels are used to determine menopause, pinpoint ovulation and monitor endocrine therapy.

PRINCIPLE OF THE ASSAY

Qualisa™ LH Quantitative Test Kit is a sandwich-based enzyme-linked immunosorbent assay. The test employs anti-LH antibody for solid phase immobilization and mouse monoclonal anti-LH antibody in the antibody-enzyme (horseradish peroxidase) conjugate solution. The test sample is allowed to react simultaneously with the two antibodies, resulting in the LH molecules being sandwiched between the solid phase and enzyme-linked antibodies. After incubation the wells are washed and labeled antibodies are detected by adding solution of TMB. The reaction is stopped after specified time with stop solution and absorbance is determined for each well using an ELISA reader. The concentration of LH is directly proportional to the color intensity of the test sample.

MATERIALS AND COMPONENTS

Materials provided with the test kits:

- Coated Microwells: Microwells coated with anti- LH antibody.
- LH Enzyme Conjugate. Ready to use.
- TMB Substrate. Ready to use
- Stop Solution. Ready to use
- LH Standard set of 6 standards labeled as A to F in lyophilized form. Ready to use. For standard Concentrations refer vial label.
- Wash Buffer Concentrate (20X).

Materials required but not provided

- Precision pipettes: 10-100µl, 20-200µl, 100-1000µl
- Disposable pipette tips
- Distilled water
- Disposable Gloves
- ELISA reader
- ELISA washer

STORAGE AND STABILITY

1. **Qualisa™ LH** kit is stable at 2-8°C upto expiry date printed on the label.
2. Coated microwells should be used within one month upon opening the pouch provided that once opened, the pouch must be resealed to protect from moisture. If the colour of the dessicant has changed from blue to white at the time of opening the pouch, another coated microwells pouch should be used.
3. Diluted Wash Buffer is stable upto one week when stored at 2-8°C.

SPECIMEN COLLECTION

1. Collect Blood specimen by venipuncture according to the standard procedure.
2. Only serum should be used.
3. Avoid grossly hemolytic, lipemic or turbid samples.
4. Preferably use fresh samples. However specimens can be stored up to 48 hours at 2-8°C, for short duration.
5. For longer storage, specimens can be frozen at -20°C. Thawed samples must be mixed prior to testing.
6. Do not heat inactivate before use.
7. Specimen containing precipitate or particulate matter should be clarified by centrifugation prior to use.
8. Specimen should be free from particulate matter and microbial contamination.

PRECAUTIONS


1. Bring all reagents and specimen to room temperature before use.
2. Do not pipette any material by mouth.
3. Do not eat, drink or smoke in the area where testing is done.
4. Use protective clothing and wear gloves when handling samples.
5. Use absorbent sheet to cover the working area.
6. Immediately clean up any spills with sodium hypochlorite.
7. All specimens and standards should be considered potentially infectious and discarded appropriately.
8. Neutralize acid containing waste before adding hypochlorite.
9. Do not use kit after the expiry date.
10. Do not mix components of one kit with another.
11. Always use new tip for each specimen and reagent.
12. Do not allow liquid from one well to mix with other wells.
13. Do not let the strips dry in between the steps.

REAGENT PREPARATION

1. All reagents should be brought to room temperature (18-25°C) and mixed by gently inverting or swirling prior to use. Do not induce foaming.
2. Dilute Wash Buffer 20 times (for example add 5ml concentrated buffer to 95 ml distilled or deionized water). Mix well before use.
3. Since the reference standards are lyophilized, reconstitute each standard with 0.5 ml distilled water. Allow the reconstituted material to stand for at least 20 minutes. Reconstituted standards should be sealed and stored at 2-8°C.

TEST PROCEDURE

1. Secure the desired number of coated wells in the holder. Dispense **25 µl** of Standards and Serums into the appropriate wells.
2. Dispense 50 µl of Enzyme Conjugate reagent into each well. Incubate at room temperature (18-25°C), for **45 minutes**.
3. Remove the incubation mixture by emptying the plate content into a waste container. Rinse and empty the microtiter plate 5 times with Wash Buffer (1X). Strike the microtiter plate sharply onto the absorbent paper or paper towels to remove all residual water droplets.
4. Dispense **100 µl** of TMB Substrate into each well. Incubate at room temperature (18-25°C) in the dark, for **20 minutes**.
5. Stop the reaction by adding **100 µl** of Stop Solution to each well. Gently mix for 10 seconds until the blue color completely changes to yellow.
6. Read the optical density at 450/630 nm with a microtiter plate reader within **15 minutes**.



CALCULATION OF RESULTS


Construct a standard curve by plotting the absorbance obtained from each reference standards against its concentrations in mIU/ml on the graph paper, with absorbance values on the vertical or Y axis and concentrations on the horizontal or X axis. Use the absorbance values for each specimen to determine the corresponding concentration of LH in mIU/ml from the standard curve. Any diluted specimens must be corrected by the appropriate dilution factor.

Example of Standard curve

Results of a typical standard run with optical density reading at 450nm (ref 600-700nm) shown in the Y axis against LH concentrations shown in the X axis.

Suggest: Use 4-Parameter Standard curve to calculate sample values.

LH Values (mIU/ml)	Absorbance
A	0.010
B	0.084
C	0.300
D	0.620
E	1.231
F	1.820



This standard curve is for the purpose of illustration only and should not be used to calculate samples. Each user should obtain his or her own standard curve and data.

Expected values and sensitivity

Each laboratory must establish its own normal ranges based on patient population. The results provided below are based on randomly selected clinical laboratory samples.

Adult male	0~25 mIU/ml
Female Follicular	0~40 mIU/ml
Mid-cycle	40~150 mIU/ml
Luteal	0~30 mIU/ml
Post Menopausal	20~200 mIU/ml
Pre pubertal, female	0~9 mIU/ml
Pre pubertal, male	0~17 mIU/ml

The minimum detectable concentration of human Lutenizing Hormone by this assay is estimated to be 2 mIU/ml.

PERFORMANCE CHARACTERISTICS

A) Internal Evaluation:

1. Accuracy: In an internal study **Qualisa™ LH** was evaluated against commercially available licensed kit with 90 random clinical samples, & **Qualisa™ LH** has demonstrated >98% clinical correlation with the commercially available licensed kit.