

FIAcheck™

NT-proBNP

Fluorescence Immunoassay for Quantitative detection of NT-proBNP in Human Serum/Plasma

DEVICE

FOR IN VITRO DIAGNOSTIC USE ONLY

Store at 4°C to 30°C

INTENDED USE

FIAcheck™ NT-proBNP Fluorescence Immunoassay is intended for the in-vitro quantitative measurement of NT-proBNP in human serum/plasma.

INTRODUCTION

N-terminal pro-brain natriuretic peptide (NT-proBNP) is produced predominantly by the cardiac ventricular myocytes and is released in response to myocardial stress and filling pressure and is involved in maintaining intravascular volume homeostasis. After stimulation of heart muscle cells, the natriuretic peptides are produced as prohormones (proBNP) and this is cleaved into two fragments which are secreted into the blood stream as the 32 amino acids active BNP and the N-terminal fragment of 76 amino acids designated as NT-proBNP. NT-proBNP immunoassays are widely used and are now considered to be a useful marker and have a high degree of diagnostic accuracy in clinical practice and cardiovascular research as a diagnostic tool for the occurrence and severity of heart failure (HF). Therefore NT-proBNP measurements in human blood are helpful not only for the cardiac disease diagnosis but also for evaluation of patients with suspected HF and assessment of severity of the disease.

PRINCIPLE OF THE TEST

FIAcheck™ NT-proBNP is based on principle of agglutination of antibodies/anti-sera with respective antigen in immuno-chromatographic format using fluorophores as signal generators. The **FIAcheck™ NT-proBNP** test device is coated with immobilized mouse anti-human NT-proBNP monoclonal antibody 1 on the test line, goat anti-chicken IgY in control line and a mixture of mouse anti-human NT-proBNP monoclonal antibody 2 and Chicken IgY labeled with fluorescent microspheres on the binding pad.

NT-proBNP in sample binds to the mouse anti-human NT-proBNP monoclonal antibody 2 labeled with fluorescent microspheres in the binding pad.

The fluorescent labeled Ag-Ab complex moves forward due to capillary action and is captured by the immobilized mouse anti-human NT-proBNP monoclonal antibody 1 forming a double antibody sandwich and produces the test line. Chicken IgY labeled with fluorescent microspheres binds with goat anti-chicken IgY to produce the control line. When the **FIAcheck™** test device is inserted in the **FIAcheck™** analyzer, it scans both the test line and control line. The ratio of the two fluorescence values is used to calculate the concentration of the analyte present in the sample.

MATERIALS AND COMPONENTS

Materials provided with the test kits:

- **FIAcheck™ NT-proBNP** test device in a sealed pouch with desiccant.
- QR Code card for calibration.

Materials required but not provided

- Precision pipettes: 100µl
- Disposable pipette tips
- Disposable Gloves
- **FIAcheck™** Analyzer (Time Resolved Fluorescence Immunoassay Analyzer)
- Digital Thermometer
- Stopwatch













STORAGE AND STABILITY

1. **FIAcheck™ NT-proBNP** kit is stable at 4°C to 30°C upto expiry date printed on the label. DO NOT FREEZE.
2. **FIAcheck™** Test device should be used within 30 minutes once the foil pouch is opened.
3. If the colour of the desiccant has changed from blue to pink or colourless at the time of opening the pouch, kindly discard the device and use another device.

SAMPLE COLLECTION

1. Only human Serum/Plasma samples should be used. Other bodily fluids and samples may not give accurate result.

SYMBOL KEYS

| | | | | | | | |
|---|-----------------------------------|---|------------------------------------|---|---------------------------|---|--------------|
|  | Temperature Limitation |  | Consult Instructions for use |  | Date of Manufacture |  | Do not reuse |
|  | Manufacturer |  | In vitro Diagnostic Medical Device |  | This side up |  | Use by |
|  | Contains sufficient for <n> tests |  | Catalogue Number |  | Batch Number / Lot Number |  | Device |



Manufactured by:

Zephyr Biomedicals

A Division of Tulip Diagnostics (P) Ltd.

M 46-47, Phase III B, Verna Industrial Estate, Verna, Goa - 403 722, INDIA.

Regd. Office: Gitanjali, Tulip Block, Dr. Antonio Do Rego Bagh, Alto Santacruz,

Bambolim Complex P.O., Goa - 403 202, INDIA.

Website: www.tulipgroup.com Email: sales@tulipgroup.com


2. Plasma can be anti-coagulated with Heparin and Sodium citrate or Tri sodium citrate under aseptic conditions.
3. The test should be performed within 4 hours after the sample collection at room temperature.
4. Avoid grossly hemolytic, lipemic or turbid samples.
5. Preferably use fresh samples. However, samples can be stored for 3 days at 2°C to 8°C, and if more than 3 days, they should be stored at -20°C.
6. The sample should be recovered to room temperature (18°C to 30°C) before testing. Avoid repeated freezing and thawing of samples as it can affect test values.
7. Samples containing precipitate or particulate matter should be clarified by centrifugation prior to use.
8. Samples should be free from particulate matter and microbial contamination.

PRECAUTIONS

1. Only for In vitro diagnostic use.
2. After the test device is removed from the sealed pouch, it should be used immediately or within 30 minutes of opening the pouch.
3. Do not reuse the tested **FIcheck™** device.
4. Do not use damaged **FIcheck™** test device or pouch.
5. All samples should be considered potentially infectious and discarded appropriately as per Standard Bio-Safety guidelines.
6. Do not use kit after the expiry date.
7. Do not mix components of one kit with another.
8. Always use new tip for each sample.
9. Scan QR code card specific to the lot you are using.
10. Ambient temperature of testing environment directly impacts the accuracy of results. Ideal working temperature is 18°C to 30°C.
11. The **FIcheck™** test device should be read immediately after the specified reaction time. Delay in reading might affect the accuracy of results.
12. The **FIcheck™** test device should be used only in conjunction with **FIcheck™** analyzer for accurate and reliable results.

TEST PROCEDURE

1. To calibrate the **FIcheck™ NT-proBNP** kit, scan the QR code card provided with the kit.
2. Remove **FIcheck™ NT-proBNP** test device from sealed pouch and place it horizontally on a clean table, label the device with sample identity.
3. Dispense **100 µl** of sample onto the sample port of **FIcheck™ NT-proBNP** test device.
4. Incubate at room temperature (18°C to 30°C) for **15 minutes**.
5. After 15 minutes, insert the test device immediately into the **FIcheck™** analyzer and read results.



Expected Range

Cut-Off Value:

<75 years old: 125pg/ml
 ≥75 years old: 531pg/ml

The cut-off value for NT-proBNP was determined by testing samples from 300 apparently healthy individuals. The 95th percentile of the concentration for NT-proBNP is 125pg/mL, and the 97.5th percentile of the concentration for NT-proBNP is 531pg/mL.

It is recommended that each laboratory establish its own reference range for the population it serves.

PERFORMANCE CHARACTERISTICS

1. Measuring Range: 50-25000pg/mL.
2. Lower Detection Limit: ≤50pg/mL.
3. Upper Detection Limit: ≥25000pg/mL.
4. Accuracy: Based on comparison experiments, the relative standard deviation of ≤15%, and the correlation coefficient of $r \geq 0.990$ was observed.
5. Within-Run Precision: ≤15%.
6. Between-Run Precision: ≤15%.
7. Hook Test: No hook effect with high concentration sample. Hook test was conducted with reference material exceeding the upper limit of measuring range, and the detection result was greater than the upper limit of detection.
8. In an internal Study, **FIcheck™ NT-proBNP** was evaluated against commercially available licensed kit with 100 random clinical samples and **FIcheck™ NT-proBNP** has demonstrated 100% clinical correlation with the commercially available licensed kit.

| NT-proBNP Levels | No. of samples | FIcheck™ NT-proBNP | EIA NT-proBNP |
|------------------|----------------|---------------------------|---------------|
| Normal | 80 | 80 | 80 |
| Low | 8 | 8 | 8 |
| High | 12 | 12 | 12 |

9. In an external Study, **FIcheck™ NT-proBNP** has been evaluated by a NABL accredited lab against their reference method. In this evaluation **FIcheck™ NT-proBNP** has demonstrated 100% correlation with the reference method.

*Data file: Zephyr Biomedicals (A Division of Tulip Diagnostics (P) Ltd).

LIMITATIONS OF THE ASSAY

1. As with all diagnostic tests, a definite clinical diagnosis should not be based on the results of a single test, but should only be made by the physician after all clinical and laboratory findings have been evaluated.
2. This kit is only for human serum and plasma.

WARRANTY

This product is designed to perform as described on the label and package insert. The manufacturer disclaims any implied warranty of use and sale for any other purpose. In the event of performance changes or product malfunction, please contact manufacturer.

BIBLIOGRAPHY

1. Jame LI, Carlos AC, Saif A, et al. The N-terminal pro-BNP investigation of dyspnea in the emergency department (PRIDE) study. Am J Cardiol, 2005, 95:948-954.
2. Charlotte K, Bjorn G, Lars K, et al. N-Terminal pro-B-typenatriuretic peptide and long-term mortality in stable coronary heart disease. N Engl J Med, 2005, 352:666-675.
3. Paulo B, Ana A, Joana P, et al. N-terminal-Pro-brain natriuretic peptide predicts outcome after hospital discharge in heart failure patients. Circulations, 2004, 110:2168-2174.
4. Lene SN, Jens S, Niels AK, et al. N-terminal pro-brain natriuretic peptide for discriminating between cardiac and non-cardiac dyspnoea. Eur J Heart Failure, 2004, 6:63-70.