

FIAcheck™

PCT

Fluorescence Immunoassay for Quantitative determination of Procalcitonin (PCT) in Human Serum and Plasma

DEVICE

FOR IN VITRO DIAGNOSTIC USE ONLY
Store at 4°C to 30°C

INTENDED USE

FIAcheck™ **PCT** Fluorescence Immunoassay is intended for the in-vitro quantitative measurement of Procalcitonin (PCT) in human serum and plasma.

INTRODUCTION

Procalcitonin (PCT) is considered a promising biomarker for the early detection of systemic bacterial infection resulting into sepsis and septic shock. PCT is a 116-amino acid peptide precursor of the hormone calcitonin and has a molecular mass of 13kDa. PCT levels in human serum have been proved to have a direct correlation with bacterial infection and sepsis but do not elevate in the case of viral infections. Other inflammatory markers like C-reactive protein (CRP), lack the specificity required to diagnose bacterial versus non-bacterial infections. PCT level in the blood of healthy individual is as low as 0.05ng/ml, but increases rapidly in response to a systemic bacterial infection to exceed 0.5ng/ml. The half- life of circulating PCT is short and may last for only 25-30 hours. PCT levels in the human serum have been shown to decrease following the use of appropriate antibiotic therapies. Therefore, the prognostic value of PCT has shown clinical significance by providing physicians with a positive correlation between disease severity and elevated serum PCT levels. PCT has been used to efficiently differentiate patients with Systemic Inflammatory Response Syndrome (SIRS) from those with sepsis, when compared with other biomarkers such as IL-2, IL-6, IL-8, CRP and TNF- α . Hence, PCT assays can be effectively used to accurately determine if a bacterial species as the cause of a patients systemic inflammatory reaction.

PRINCIPLE OF THE TEST

FIAcheck™ **PCT** is based on principle of agglutination of antibodies/anti-sera with respective antigen in immunochromatographic format using fluorophores as signal generators. The **FI**Acheck™ **PCT** test device is coated with immobilized PCT mouse monoclonal antibody 1 on the test line, goat anti chicken IgY in control line and a mixture of PCT mouse monoclonal antibody 2 and Chicken IgY labeled with fluorescent microspheres on the binding pad. PCT in sample binds to the PCT mouse monoclonal antibody 2 labeled with fluorescent microspheres in the binding pad. The fluorescent labeled Ag-Ab complex moves forward due to capillary action and is captured by the immobilized PCT mouse monoclonal antibody 1 forming a double antibody sandwich and produces the test line. Chicken IgY labelled with fluorescent microspheres binds with goat anti chicken IgY to produce the control line. When the **FI**Acheck™ test device is inserted in the **FI**Acheck™ analyzer, it scans both the test line and control line. The ratio of the two fluorescence values is used to calculate the concentration of the analyte present in the sample.

MATERIALS AND COMPONENTS

Materials provided with the test kits:

- **FI**Acheck™ **PCT** test device in a sealed pouch with desiccant.
- QR Code card for calibration.













Materials required but not provided

- Precision pipettes: 100 μ l
- Disposable pipette tips
- Disposable Gloves
- **FI**Acheck™ Analyzer (Time Resolved Fluorescence Immunoassay Analyzer)
- Digital Thermometer
- Stopwatch

STORAGE AND STABILITY

1. **FI**Acheck™ **PCT** kit is stable at 4°C to 30°C upto expiry date printed on the label. DO NOT FREEZE.
2. **FI**Acheck™ Test device should be used within 30 minutes once the foil pouch is opened.
3. If the colour of the desiccant has changed from blue to pink or colourless at the time of opening the pouch, kindly discard the device and use another device.

SYMBOL KEYS

 Temperature Limitation	 Consult Instructions for use	 Date of Manufacture	 Do not reuse
 Manufacturer	 IVD In vitro Diagnostic Medical Device	 This side up	 Use by
 Contains sufficient for <n> tests	 REF Catalogue Number	 LOT Batch Number / Lot Number	 DEVICE Device



Manufactured by:

Zephyr Biomedicals

A Division of Tulip Diagnostics (P) Ltd.

M 46-47, Phase III B, Verna Industrial Estate, Verna, Goa - 403 722, INDIA.

Regd. Office: Gitanjali, Tulip Block, Dr. Antonio Do Rego Bagh, Alto Santacruz, Bambolim Complex P.O., Goa - 403 202, INDIA.

SAMPLE COLLECTION


1. Only human Serum and Plasma samples should be used. Other bodily fluids and samples may not give accurate result.
2. Plasma can be anti-coagulated with Heparin and Sodium citrate or Tri sodium citrate under aseptic conditions.
3. The test should be performed within 4 hours after the sample collection at room temperature.
4. Avoid grossly hemolytic, lipemic or turbid samples.
5. Preferably use fresh samples. However, samples can be stored for 3 days at 2°C to 8°C, and if more than 3 days, they should be stored at -20°C.
6. The sample should be recovered to room temperature (18°C to 30°C) before testing. Avoid repeated freezing and thawing of samples as it can affect the test values.
7. Samples containing precipitate or particulate matter should be clarified by centrifugation prior to use.
8. Samples should be free from particulate matter and microbial contamination.

PRECAUTIONS

1. Only for In vitro diagnostic use.
2. After the test device is removed from the sealed pouch, it should be used immediately or within 30 minutes of opening the pouch.
3. Do not reuse the tested **FIcheck™** device. Do not use sample dilution vial for more than one sample.
4. Do not use damaged **FIcheck™** test device or pouch.
5. All samples should be considered potentially infectious and discarded appropriately as per Standard Bio-Safety guidelines.
6. Do not use kit after the expiry date.
7. Do not mix components of one kit with another.
8. Always use new tip for each sample and reagent.
9. Scan QR code card specific to the lot you are using.
10. Ambient temperature of testing environment directly impacts the accuracy of results. Ideal working temperature is 18°C to 30°C.
11. The **FIcheck™** test device should be read immediately after the specified reaction time. Delay in reading might affect the accuracy of results.
12. The **FIcheck™** test device should be used only in conjunction with **FIcheck™** analyzer for accurate and reliable results.

TEST PROCEDURE

1. To calibrate the **FIcheck™ PCT** kit, scan the QR code card provided with the kit.
2. Remove **FIcheck™ PCT** test device from sealed pouch and place it horizontally on a clean table, label the device with sample identity.
3. Dispense 100 µl of sample onto the sample port of **FIcheck™ PCT** test device.
4. Incubate at room temperature (18°C to 30°C) for 10 minutes.
5. After 10 minutes, insert the test device immediately into the **FIcheck™** analyzer and read results.



Expected Range

Cut-Off Value: 0.5ng/mL

PCT concentration is determined using samples obtained from 200 apparently healthy individuals. It is recommended that each laboratory establish its own reference range for the population it serves.

PERFORMANCE CHARACTERISTICS

1. Measuring Range: 0.1-50ng/mL.
2. Lower Detection Limit: ≤ 0.1 ng/mL.
3. Upper Detection Limit: ≥ 50 ng/mL.
4. Accuracy: Based on comparison experiments, the relative standard deviation of $\leq 15\%$, and the correlation coefficient of $r \geq 0.990$ was observed.
5. Within-Run Precision: $\leq 15\%$.
6. Between-Run Precision: $\leq 15\%$.
7. Hook Test: No hook effect with high concentration sample. Hook test was conducted with reference material exceeding the upper limit of measuring range, and the detection result was greater than the upper limit of detection.
8. In an internal Study, **FIcheck™ PCT** was evaluated against commercially available licensed kit with 100 random clinical samples and **FIcheck™ PCT** has demonstrated 100% clinical correlation with the commercially available licensed kit.

PCT Levels	No. of samples	FIcheck™ PCT	EIA PCT
Normal	80	80	80
Low	8	8	8
High	12	12	12

9. In an external Study, **FIcheck™ PCT** has been evaluated by a NABL accredited lab against their reference method. In this evaluation **FIcheck™ PCT** has demonstrated 100% correlation with the reference method.

*Data file: Zephyr Biomedicals (A Division of Tulip Diagnostics (P) Ltd).

LIMITATIONS OF THE ASSAY

1. As with all diagnostic tests, a definite clinical diagnosis should not be based on the results of a single test but should only be made by the physician after all clinical and laboratory findings have been evaluated.
2. This kit is only for human serum and plasma.

WARRANTY

This product is designed to perform as described on the label and package insert. The manufacturer disclaims any implied warranty of use and sale for any other purpose. In the event of performance changes or product malfunction, please contact manufacturer.

BIBLIOGRAPHY

1. Lorrot M, Morhn F, Coste J, et al. Procalcitonin in pediatric emergencies: comparison with C-reactive protein. Interleukin and interferon alpha in the differentiation between bacteria and viral infections. Presse Medicale, 2000, 29(3): 128-34.
2. Guven H, Aldntop L, Baydin A, et al. Diagnostic value of procalcitonin levels as an early indicator of sepsis. Am J EmergMed, 2002, 20(3): 202-206.
3. Balei C, Sungurtekin H, Gttrses E, et al. Usefulness of procalcitonin for diagnosis of sepsis in the intensive care unit. Crit Care. 2003, 7(1): 85-90.
4. Yukioka H, Yoshida G, Kurita S. Plasma procalcitonin in sepsis and organ failure. Ann Acad Med Singapore, 2001, 30(5): 528-531.