

# **FI**Acheck™

## **PRL**

**Fluorescence Immunoassay for Quantitative determination of Prolactin (PRL) in Human Serum**

**DEVICE**

FOR IN VITRO DIAGNOSTIC USE ONLY

Store at 4°C to 30°C

### INTENDED USE

**FI**Acheck™ PRL Fluorescence Immunoassay is intended for the in-vitro quantitative measurement of Prolactin (PRL) in human serum.

### INTRODUCTION

Human prolactin (lactogenic hormone) is secreted from the anterior pituitary gland in both men and women. Women normally have slightly higher basal prolactin levels than men; apparently, there is an estrogen-related rise at puberty and a corresponding decrease at menopause. The primary functions of prolactin are to initiate breast development and to maintain lactation. Prolactin also suppresses gonadal function. The determination of prolactin concentration is helpful in diagnosing hypothalamic-pituitary disorders. High prolactin levels are commonly associated with galactorrhea and amenorrhoea. Prolactin concentrations have been shown to be increased by estrogens, thyrotropin-releasing hormone (TRH) and several drugs affecting dopaminergic mechanism. Prolactin levels are elevated in renal disease and hypothyroidism and in some situations of stress, exercise, and hypoglycemia.

Prolactin concentrations may also be increased by drugs such as chlorpromazine and reserpine and may be lowered by bromocryptine and L-dopa.

### PRINCIPLE OF THE TEST

**FI**Acheck™ PRL is based on principle of agglutination of antibodies/anti-sera with respective antigen in immunochromatographic format using fluorophores as signal generators. The **FI**Acheck™ PRL test device is coated with immobilized PRL mouse monoclonal antibody 1 on the test line, sheep anti chicken IgY in control line and a mixture of PRL mouse monoclonal antibody 2 and Chicken IgY labeled with fluorescent microspheres on the binding pad. PRL in sample binds to the PRL mouse monoclonal antibody 2 labeled with fluorescent microspheres in the binding pad. The fluorescent labeled Ag-Ab complex moves forward due to capillary action and is captured by the immobilized PRL mouse monoclonal antibody 1 forming a double antibody sandwich and produces the test line. Chicken IgY labelled with fluorescent microspheres binds with sheep anti chicken IgY to produce the control line. When the **FI**Acheck™ test device is inserted in the **FI**Acheck™ analyzer, it scans both the test line and control line. The ratio of the two fluorescence values is used to calculate the concentration of the analyte present in the sample.

### MATERIALS AND COMPONENTS

#### Materials provided with the test kits:

- **FI**Acheck™ PRL test device in a sealed pouch with desiccant.
- QR Code card for calibration.

#### Materials required but not provided

- Precision pipettes: 100µl
- Disposable pipette tips
- Disposable Gloves
- **FI**Acheck™ Analyzer (Time Resolved Fluorescence Immunoassay Analyzer)
- Digital Thermometer
- Stopwatch













### STORAGE AND STABILITY

1. **FI**Acheck™ PRL kit is stable at 4°C to 30°C upto expiry date printed on the label. DO NOT FREEZE.
2. **FI**Acheck™ Test device should be used within 30 minutes once the foil pouch is opened.
3. If the colour of the desiccant has changed from blue to pink or colourless at the time of opening the pouch, kindly discard the device and use another device.

### SAMPLE COLLECTION

1. Only Human Serum sample should be used. Other bodily fluids and samples may not give accurate result.
2. The test should be performed within 4 hours after the sample collection at room temperature.
3. Avoid grossly hemolytic, lipemic or turbid samples.

### SYMBOL KEYS

 Temperature Limitation	 Consult Instructions for use	 Date of Manufacture	 Do not reuse
 Manufacturer	 IVD In vitro Diagnostic Medical Device	 This side up	 Use by
 Contains sufficient for <n> tests	 REF Catalogue Number	 LOT Batch Number / Lot Number	 DEVICE Device



Manufactured by:

**Zephyr Biomedicals**

A Division of Tulip Diagnostics (P) Ltd.

M 46-47, Phase III B, Verna Industrial Estate, Verna, Goa - 403 722, INDIA.

Regd. Office: Gitanjali, Tulip Block, Dr. Antonio Do Rego Bagh, Alto Santacruz, Bambolim Complex P.O., Goa - 403 202, INDIA.


4. Preferably use fresh samples. However, samples can be stored for 3 days at 2°C to 8°C, and if more than 3 days, they should be stored at -20°C.
5. The sample should be recovered to room temperature (18°C to 30°C) before testing. Avoid repeated freezing and thawing of samples as it can affect the test values.
6. Samples containing precipitate or particulate matter should be clarified by centrifugation prior to use.
7. Samples should be free from particulate matter and microbial contamination.

#### PRECAUTIONS

1. Only for In vitro diagnostic use.
2. After the test device is removed from the sealed pouch, it should be used immediately or within 30 minutes of opening the pouch.
3. Do not reuse the tested **FIcheck™** device.
4. Do not use damaged **FIcheck™** test device or pouch.
5. All samples should be considered potentially infectious and discarded appropriately as per Standard Bio-Safety guidelines.
6. Do not use kit after the expiry date.
7. Do not mix components of one kit with another.
8. Always use new tip for each sample.
9. Scan QR code card specific to the lot you are using.
10. Ambient temperature of testing environment directly impacts the accuracy of results. Ideal working temperature is 18°C to 30°C.
11. The **FIcheck™** test device should be read immediately after the specified reaction time. Delay in reading might affect the accuracy of results.
12. The **FIcheck™** test device should be used only in conjunction with **FIcheck™** analyzer for accurate and reliable results.

#### TEST PROCEDURE

1. To calibrate the **FIcheck™ PRL** kit, scan the QR code card provided with the kit.
2. Remove **FIcheck™ PRL** test device from sealed pouch and place it horizontally on a clean table, label the device with sample identity.
3. Dispense 100 µl of sample onto the sample port of **FIcheck™ PRL** test device.
4. Incubate at room temperature (18°C to 30°C) for 15 minutes.
5. After 15 minutes, insert the test device immediately into the **FIcheck™** analyzer and read results.



#### Expected Range

The following reference interval was obtained after statistical analysis of the confidence interval for the tests of the content of PRL in serum samples of healthy people:

Mature male: 86.30-425.72µIU/mL

Mature female: 72.55-600.4µIU/mL

It is recommended that each laboratory establish its own reference range for the population it serves.

#### PERFORMANCE CHARACTERISTICS

1. Measuring Range: 25-4000µIU/mL.
2. Lower Detection Limit: ≤20µIU/mL.
3. Upper Detection Limit: ≥4000µIU/mL.
4. Accuracy: Based on comparison experiments, the relative standard deviation of ≤15%, and the correlation coefficient of  $r \geq 0.990$  was observed.
5. Within-Run Precision: ≤15%.
6. Between-Run Precision: ≤15%.
7. Hook Test: No hook effect with high concentration sample. Hook test was conducted with reference material exceeding the upper limit of measuring range, and the detection result was greater than the upper limit of detection.
8. In an internal Study, **FIcheck™ PRL** was evaluated against commercially available licensed kit with 100 random clinical samples and **FIcheck™ PRL** has demonstrated 100% clinical correlation with the commercially available licensed kit.

PRL Levels	No. of samples	<b>FIcheck™ PRL</b>	EIA PRL
Normal	80	80	80
Low	8	8	8
High	12	12	12

9. In an external Study, **FIcheck™ PRL** has been evaluated by a NABL accredited lab against their reference method. In this evaluation **FIcheck™ PRL** has demonstrated 100% correlation with the reference method.

\*Data file: Zephyr Biomedicals (A Division of Tulip Diagnostics (P) Ltd).

#### LIMITATIONS OF THE ASSAY

1. As with all diagnostic tests, a definite clinical diagnosis should not be based on the results of a single test, but should only be made by the physician after all clinical and laboratory findings have been evaluated.
2. This kit is only for human serum.

#### WARRANTY

This product is designed to perform as described on the label and package insert. The manufacturer disclaims any implied warranty of use and sale for any other purpose. In the event of performance changes or product malfunction, please contact manufacturer.

#### BIBLIOGRAPHY

1. Yang Yue. Clinical significance of measuring six sex hormones in women[J]. Chinese Journal of Clinicians, 2003, 31(4):50-51.
2. Scott MG, Ladenson JH, Green ED, Gast MJ. Hormonal evaluation of female infertility and reproductive disorders. Clin Chem. 1989, 35: 620– 630.
3. Huerta R, Malacara JM, Fajardo ME, Nava LE, Bocanegra A, Sanchez J. High-Frequency FSH and PRL Pulses in Obese Menopausal Women. Endocrine, 1997, 7 (3): 281-286.
4. Pincus SM, Veldhuis JD, Mulligan H, Iranmanish A, Evans WS. Effects of age on the irregularity of PRL and FSH serum concentrations in women and men. Am J Physiol 273 (Endocrinol Metab 36), 1997, E989-E995.