



parascreen®

| REF | 503030010 | 503030025 | 503030050 | 503030100 |
|-----|-----------|-----------|-----------|-----------|
|-----|-----------|-----------|-----------|-----------|

## RAPID TEST FOR MALARIA

Pan/Pf

DEVICE

### INTENDED USE

**parascreen®** is a rapid, qualitative, two site sandwich immunoassay utilizing capillary and venous whole blood specimens of symptomatic patients for the detection of *P.falciparum* specific histidine rich protein-2 (Pf. HRP-2) and Plasmodium Lactate Dehydrogenase (pLDH) antigens produced by *P.falciparum*, *P.vivax*, *P.ovale* and *P.malariae* species. It is used in the diagnosis of malaria for differentiation of *P.falciparum* and other malaria species.

**parascreen®** is intended to be used by trained healthcare or laboratory professionals or other health care workers who have received appropriate training. This product can be used by trained lay providers operating at point-of-care in resource-limited settings. This product is not intended for self-testing and it is not for blood donor screening. The test is not automated; it needs to be performed and interpreted manually by the user.

### SUMMARY

Four species of the Plasmodium parasites are responsible for malaria infections in human viz. *P. falciparum*, *P.vivax*, *P.ovale* and *P.malariae*. Of these, *P. falciparum* and *P.vivax* are the most prevalent. Early detection and differentiation of malaria is of utmost importance due to incidence of cerebral malaria and drug resistance associated with falciparum malaria and due to the morbidity associated with the other malarial forms.

**parascreen®** detects the presence of Pan malaria specific pLDH released from parasitised blood cells, for the detection of all malarial parasites. Whereas, for the detection of *P. falciparum* malaria, **parascreen®** utilises the detection of *P. falciparum* specific histidine rich protein-2 (Pf. HRP-2) which is a water soluble protein that is released from parasitised red blood cells of infected individuals. In the absence of *P.falciparum* specific Pf. HRP-2, the presence of Pan malaria specific band points to the presence of other malarial species such as *P.vivax*, *P.ovale* or *P.malariae*. Speciation is done and results inferred in the context of prevalence rates of the malarial species prevalent in the particular region.

### PRINCIPLE

**parascreen®** utilizes the principle of agglutination of antibodies/ antisera with respective antigen in immuno-chromatography format along with use of nano gold particles as agglutination revealing agent. As the test specimen flows through the membrane assembly of the device after addition of the clearing buffer, the colored colloidal gold conjugates of the Agglutinating sera for HRP-2 / Agglutinating sera for Pan malaria specific pLDH - colloidal gold conjugate complexes the proteins in the lysed specimen. This complex moves further on the membrane to the test region where it is immobilised by the Agglutinating sera for HRP-2 / Agglutinating sera for Pan malaria specific pLDH coated on the membrane leading to formation of pink-purple colored band/s which confirms a positive test result. Absence of this colored band/s in the test region indicates a negative test result.

The unreacted conjugate along with the rabbit globulin-colloidal gold conjugate and unbound complex if any, move further on the membrane and are subsequently immobilised by Agglutinating sera for rabbit globulin coated on the membrane at the control region, forming a pink-purple band. The control band formation is based on the 'Rabbit / Agglutinating sera for Rabbit globulin' system. Since it is completely independent of the analyte detection system, it facilitates formation of consistent control band signal independent of the analyte concentration. This control band serves to validate the test performance.

### REAGENTS AND MATERIALS SUPPLIED

**parascreen®** kit contains:

- Individual pouches, each containing:
  - DEVICE** Test Device: Membrane assembly pre-dispensed with Agglutinating sera for HRP-2 - colloidal gold conjugate, Agglutinating sera for Pan malaria specific pLDH - colloidal gold conjugate, rabbit globulin-colloidal gold conjugate, Agglutinating sera for HRP-2, Agglutinating sera for Pan malaria specific pLDH and Agglutinating sera for rabbit globulin at the respective regions.
  - Desiccant pouch.
  - PIPETTE** Disposable Plastic Specimen Applicator.
- BUF** Clearing Buffer in a dropper bottle.
- Instructions for use.
- Pictorial instructions for use.
- Alcohol swabs – 70% Isopropyl alcohol
- Sterile lancets.

| Product codes                  | REF | 503030010  | 503030025  | 503030050  | 503030100  |
|--------------------------------|-----|------------|------------|------------|------------|
| Pouch sealed tests             | ▽   | 10         | 25         | 50         | 100        |
| Clearing buffer bottles        |     | 01 x 3.0ml | 01 x 4.0ml | 02 x 4.0ml | 04 x 4.0ml |
| Alcohol swabs                  |     | 10         | 25         | 50         | 100        |
| Sterile lancets                |     | 10         | 25         | 50         | 100        |
| Instructions for use           |     | 01         | 01         | 01         | 01         |
| Pictorial instructions for use |     | 01         | 01         | 01         | 01         |

### MATERIALS REQUIRED BUT NOT PROVIDED

Calibrated micropipette capable of delivering 5µl specimen accurately, disposable micropipette tips.

Permanent marker Pen/pencil, disposable gloves, timer.

Biosafety sharps container and Biohazard waste container (for potentially infectious waste).

Venipuncture blood collection kit (if whole blood is collected by venepuncture).

Additional alcohol swabs (if any included in the kit are found dry) and additional sterile lancets (if any included in the kit have the sterility seal broken).

### STORAGE AND STABILITY

The sealed pouches in the test kit & the kit components may be stored between 4°C to 40°C till the duration of the shelf life as indicated on the pouch/ carton. DO NOT FREEZE. After first opening of the clearing buffer bottle, it can be stored between 4°C to 40°C for the remaining duration of its shelf life.

### WARNINGS

Read the instructions carefully before performing the test.

For in vitro diagnostic use only. NOT FOR MEDICINAL USE. For professional use.

The test is for aiding in diagnosis of malaria infection and not for screening which requires confirmation.

Do not use beyond expiry date.

Do not use components from different lots of the product.

The device, specimen applicator, alcohol swab and blood lancet are for single use only.

Contact with the contents of desiccant pouch containing, among other substances, cobalt chloride (CAS# 7646-79-9) should be kept to a minimum. Inhalation / swallowing may cause harm.

Handle all specimens as potentially infectious.

Follow standard biosafety guidelines for handling and disposal of potentially infectious material.

Clearing buffer contains Sodium Azide(0.1%), avoid skin contact with this reagent. Azide may react with lead and copper in the plumbing system and form highly explosive metal oxides. Flush with large volumes of water to prevent azide build up in the plumbing.

### SPECIMEN COLLECTION AND PREPARATION

**For specimen collection, refer to pictorial instructions for use.**


Fresh capillary/venous wholeblood from finger prick / puncture should be used as a test specimen. However, fresh anti-coagulated venous whole blood may also be used as a test specimen. Using standard blood collection practices, collect venous whole blood into the commercially available anti-coagulant tube such as EDTA or CPDA or Heparin or Oxalate or Tri-sodium Citrate. If immediate testing is not possible then the specimen may be stored at 2°C to 8°C for upto 72 hours before testing and should be brought to room temperature (20°C to 30°C) before use on the test. Clotted, hemolysed or lipaemic whole blood specimens should not be used for performing the test.

## TESTING PROCEDURE AND INTERPRETATION OF RESULTS

1. Bring the **parascreen**<sup>®</sup> kit components to room temperature (20°C to 30°C) before testing.
  2. Open the pouch and retrieve the device, specimen applicator and desiccant pouch. Check the color of the desiccant. It should be blue, if it has turned colorless or pink, discard the device and use another device. **Once opened, the device must be used immediately.**
  3. Label the test device with patient identifier.
  4. Place the testing device on a flat horizontal surface.
  5. Tighten the cap of the clearing buffer bottle provided with the kit in the clockwise direction to pierce the buffer bottle nozzle.
  6. **Specimen application**
    - 6.1 **Venous whole blood:** Evenly mix the anti-coagulated whole blood by gentle swirling. Dip the specimen applicator into the whole blood. Ensuring that an applicator full of blood is retrieved, immediately blot the blood so collected in the specimen port 'A' (This delivers approximately 5µl of the whole blood specimen). Alternatively, 5µl of the anti-coagulated venous whole blood specimen may be delivered in the specimen port 'A' using a micropipette.
    - 6.2 **Capillary whole blood:** Touch the specimen applicator to the whole blood on the finger prick. Ensuring that an applicator full of blood is retrieved, immediately blot the blood so collected in the specimen port 'A' (Care should be taken that whole blood specimen is not clotted and transfer to the specimen port is immediate). Alternatively, 5µl of the capillary finger-prick whole blood specimen may be delivered in the specimen port 'A' using a micropipette.
- Note:** Ensure that the whole blood from the specimen applicator has been completely taken up at the specimen port 'A'.
7. Immediately dispense **two drops** of clearing buffer into buffer port 'B' holding the buffer bottle vertically and switch on the timer. To avoid contamination of clearing buffer bottle, do not touch the buffer port 'B' with the tip of clearing buffer bottle.
  8. Read the results at the end of **20 minutes** as follows :



**NEGATIVE for malaria:** Only one pink-purple band appears in the control window 'C'.





**POSITIVE for *P. falciparum*:** In addition to the control band, two pink-purple bands appear at regions 'Pf' and 'Pan' in the test window 'T'.  
OR  
In addition to the control band, a pink-purple band appear only at region 'Pf' in the test window 'T'.  
Appearance of coloured bands of any intensity (faint to dark) at 'Pf' and/or 'Pan' should be considered as positive result for *P. falciparum*.



**POSITIVE for Mixed infection (*P. falciparum* and *P. vivax* or *P. malariae* or *P. ovale*) :** In addition to the control band, two pink-purple bands appear at regions 'Pf' and 'Pan' in the test window 'T'. Appearance of coloured bands of any intensity (faint to dark) at 'Pf' and 'Pan' should be considered as positive result for Mixed infection.



**POSITIVE for Other species (non falciparum):** In addition to the control band, one pink-purple band appears only at region 'Pan' in the test window 'T'. Appearance of a coloured band of any intensity (faint to dark) at 'Pan' should be considered as positive result *P. vivax* or *P. malariae* or *P. ovale* malaria.



**INVALID RESULT:** The test should be considered invalid if no bands appear on the device. The test should also be considered invalid if only test bands (Pan and/or Pf) appear and no control band appears. Repeat the test with a new device ensuring that the test procedure has been followed accurately.

**CAUTION :** Do not read results after 30 minutes as it may show erroneous results.

## PERFORMANCE CHARACTERISTICS

### A. Analytical Performance Study

#### A1. Potentially interfering exogenous and endogenous substances :

The following Potentially interfering substances have no impact on test results of **parascreen**<sup>®</sup> :

| Type of Specimen     |                       | Sr. No. | Potential Interfering substances |               |
|----------------------|-----------------------|---------|----------------------------------|---------------|
| Endogenous substance |                       | 1       | Total Protein                    |               |
|                      |                       | 2       | Bilirubin, conjugated            |               |
|                      |                       | 3       | Cholesterol                      |               |
|                      |                       | 4       | Triglycerides                    |               |
|                      |                       | 5       | Haemoglobin                      |               |
| Common Drugs         | Antibiotic            | 1       | Amoxicillin                      |               |
|                      |                       | 2       | Ciprofloxacin                    |               |
|                      | Anti-inflammatory     | 1       | Aspirin                          |               |
|                      |                       | 2       | Ibuprofen                        |               |
| Exogenous Substance  | Anti-Malaria Drugs    |         | 1                                | Chloroquine   |
|                      |                       |         | 2                                | Doxycycline   |
|                      |                       |         | 3                                | ACT           |
|                      |                       |         | 4                                | Primaquine    |
|                      |                       |         | 5                                | Mefloquine    |
|                      |                       |         | 6                                | Sulfadoxine   |
|                      |                       |         | 7                                | Pyrimethamine |
|                      | Anti-TB Drugs         |         | 1                                | Ethambutol    |
|                      |                       |         | 2                                | Isoniazide    |
|                      |                       |         | 3                                | Rifampin      |
|                      | Anti-Retroviral Drugs |         | 1                                | Lamivudine    |
|                      |                       |         | 2                                | Efavirenz     |
|                      |                       |         | 3                                | Emtricitabine |
|                      |                       |         | 4                                | Tenofovir     |
|                      |                       |         | 5                                | Atazanavir    |

#### A2. Cross Reacting infections, disease and medical conditions:

The following 17 potential cross reacting infections/diseases/conditions did not affect the performance of **parascreen**<sup>®</sup> .

| Potential Cross reacting infections/diseases/conditions |                               |    |                           |
|---|-------------------------------|----|---------------------------|
| 1   | <i>T. cruzi</i>               | 10 | <i>Toxoplasma gondii</i>  |
| 2   | Dengue virus                  | 11 | Influenza A/B             |
| 3   | <i>Leishmania spp</i>         | 12 | Yellow fever virus        |
| 4   | <i>Brucella spp</i>           | 13 | <i>Leptospira spp</i>     |
| 5   | Measles virus (Rubeola virus) | 14 | <i>Treponema pallidum</i> |
| 6   | HAV                           | 15 | HAMA                      |
| 7   | HBV                           | 16 | ANA                       |
| 8   | HCV                           | 17 | Rheumatoid factor         |
| 9   | HIV-1/HIV-2                   |    |                           |

A3. Precision (Repeatability)

Within run, precision was determined using 10 replicates of 5 different venous whole blood specimens in 03 different lots of **parascreen**® which is summarized below:

| *Quality control Panel                           | Accuracy (%) |
|--|--------------|
| Malaria Negative                                 | 100%         |
| <i>P.falciparum</i> Positive (Moderate Positive) | 100%         |
| <i>P.falciparum</i> Positive (Weak Positive)     | 100%         |
| <i>P.vivax</i> Positive (Moderate Positive)      | 100%         |
| <i>P.vivax</i> Positive (Weak Positive)          | 100%         |

A4. Precision (Reproducibility)

Between run, precision was determined using 5 different blinded venous whole blood specimens in 3 different lots of **parascreen**® X 3 different operators X 3 different sites X 5 different days which is summarized below:

| *Quality control Panel                           | Accuracy (%) |                  |             |              |
|--|--------------|------------------|-------------|--------------|
|  | Between Day  | Between Operator | Between lot | Between site |
| Malaria Negative                                 | 100%         | 100%             | 100%        | 100%         |
| <i>P.falciparum</i> Positive (Moderate Positive) | 100%         | 100%             | 100%        | 100%         |
| <i>P.falciparum</i> Positive (Weak Positive)     | 100%         | 97.7%            | 100%        | 100%         |
| <i>P.vivax</i> Positive (Moderate Positive)      | 100%         | 100%             | 100%        | 100%         |
| <i>P.vivax</i> Positive (Weak Positive)          | 100%         | 100%             | 100%        | 100%         |

\*Quality control panel specimens have been confirmed by microscopy as malaria negative and malaria positive. Malaria positive specimens were classified as moderate or weak positive based on respective parasite counts as determined by microscopy.

A5. Analytical Sensitivity

The sensitivity of **parascreen**® for *P.falciparum* is 100 parasites/µl and for *P.vivax* is 200 parasite/µl based on microscopy results.

B. Clinical Performance study: Diagnostic Specificity and Diagnostic Sensitivity

B1. In an in-house study, a panel of 251 venous whole blood specimens whose results were earlier confirmed with microscopy were tested with **parascreen**®. The results obtained are as follows:

| Specimens           | Total no. of specimens tested | <b>parascreen</b> ® |          | Sensitivity (95% CI)     | Specificity (95% CI)     |
|---------------------|-------------------------------|---------------------|----------|--------------------------|--------------------------|
|                     |                               | Positive            | Negative |                          |                          |
| <i>P.falciparum</i> | 16                            | 16                  | 0        | 100% (79.41% to 100.00%) | -                        |
| <i>P.vivax</i>      | 25                            | 25                  | 0        | 100% (86.28% to 100.00%) | -                        |
| Malaria Negative    | 210                           | 0                   | 210      | -                        | 100% (98.26% to 100.00%) |

B2. External evaluation studies:

Table 1

| Study Site           | Total Number of Malaria Negative specimens Tested | Specimen Type                         |                                 | Number of specimens Negative by Microscopy | Number of specimens Negative in <b>parascreen</b> ® | Number of specimens falsely Positive in <b>parascreen</b> ® |
|----------------------|---|---------------------------------------|---------------------------------|--|---|---|
|                      |   | Population type                       | Mode of Collection              |  |   |   |
| Jharkhand, India     | 985   | Hospitalized Patients                 | Finger prick/ venous phlebotomy | 985  | 985   | 0   |
| Maharashtra, India.  | 1000  | Blood Donors                          | Venous whole blood              | 1000                                       | 1000  | 0   |
| Goa, India.          | 39  | Symptomatic/ Asymptomatic Individuals | Capillary Whole Blood           | 39   | 39  | 0   |
|                      |   |                                       | Venous Whole Blood              | 39   | 39  | 0   |
| Odisha, India        | 545   | Pregnant Women                        | Venous Whole Blood              | 545  | 545   | 0   |
| Odisha, India        | 497   | Neonates                              | Heel Prick                      | 497  | 497   | 0   |
| Based on above data: |   |                                       |                                 |  |   |   |
| Total Nos. tested    |   |                                       | Overall Specificity             |  | 95% Confidence Interval                             |   |
| 3105                 |   |                                       | 100%                            |  | 99.88% to 100.00%                                   |   |

Table 2

| Study Site                           | Total Number of Malaria Positive specimens Tested | Population type                       | Specimen Type                    |                                      | Number of specimens Positive by Microscopy | Number of specimens Positive in <b>parascreen</b> ® | Number of specimens falsely Negative in <b>parascreen</b> ® |
|--------------------------------------|---|---------------------------------------|----------------------------------|--------------------------------------|--|---|---|
|                                      |   |                                       | Mode of Collection               | Species Type                         |  |   |   |
| India                                | 403   | Hospitalized Patients                 | Finger prick/ venous phlebotomy  | <i>P.falciparum</i>                  | 403  | 403   | 0   |
|                                      | 312   |                                       |                                  | <i>P.vivax</i>                       | 312  | 312   | 0   |
|                                      | 26  | Symptomatic/ Asymptomatic Individuals | Capillary and Venous Whole Blood | <i>P.falciparum</i>                  | 26   | 26  | 0   |
|                                      | 29  |                                       |                                  | <i>P.vivax</i>                       | 29   | 29  | 0   |
|                                      | 06  |                                       |                                  | <i>P.falciparum</i> + <i>P.vivax</i> | 06   | 06  | 0   |
|                                      | 01  | Pregnant Women                        | Venous Whole Blood               | <i>P.falciparum</i>                  | 01   | 01  | 0   |
|                                      | 04  |                                       |                                  | <i>P.vivax</i>                       | 04   | 04  | 0   |
|                                      | 01  | Neonates                              | Heel Prick                       | <i>P.falciparum</i>                  | 01   | 01  | 0   |
|                                      | 02  |                                       |                                  | <i>P.vivax</i>                       | 02   | 02  | 0   |
|                                      | 25  | Hospitalized Patients                 | Venous whole blood               | <i>P.ovale</i>                       | 25   | 25  | 0   |
| Malawi                               | 25  |                                       |                                  | <i>P.malariae</i>                    | 25   | 25  | 0   |
| Malaysia                             | 90  | Hospitalized Patients                 | Venous whole blood               | <i>P. knowlesi</i>                   | 90   | 74  | 16  |
| Based on above data:                 |   |                                       |                                  |                                      |  |   |   |
| Plasmodium species                   |   | Total Nos. tested                     |                                  | Overall Sensitivity                  |  | 95% Confidence Interval                             |   |
| <i>P.falciparum</i>                  |   | 431                                   |                                  | 100%                                 |  | 99.15% to 100.00%                                   |   |
| <i>P.vivax</i>                       |   | 347                                   |                                  | 100%                                 |  | 98.94% to 100.00%                                   |   |
| <i>P.ovale</i>                       |   | 25                                    |                                  | 100%                                 |  | 84.56% to 100.00%                                   |   |
| <i>P.malariae</i>                    |   | 25                                    |                                  | 100%                                 |  | 84.56% to 100.00%                                   |   |
| <i>P. knowlesi</i>                   |   | 90                                    |                                  | 82.22%                               |  | 72.74% to 89.48%                                    |   |
| <i>P.falciparum</i> + <i>P.vivax</i> |   | 06                                    |                                  | 100%                                 |  | 54.07% to 100.00%                                   |   |

LIMITATIONS OF THE TEST

- 1. As with all diagnostic tests, the test result must always be correlated with clinical findings.
- 2. The results of test are to be interpreted within the epidemiological, clinical and therapeutic context.
- 3. Any modification to the above procedure and / or use of other reagents will invalidate the test procedure.
- 4. Hook effect may be observed at parasite density ≥3,00,000 parasite/µl. In such cases, repeat the test by using different dilutions of same specimen. Other clinical data (e.g symptoms, travel history, risky factors) should be used in conjunction with the test results.
- 5. Interference due to presence of heterophile antibodies in patient's specimen can lead to erroneous analyte detection in immunoassay, has been reported in various studies. **parascreen®** uses HETEROPHILIC BLOCKING REAGENT (HBR) to inhibit majority of these interferences.
- 6. Potential cross-reacting diseases such as HAT, Tick-borne Encephalitis and those caused by *Schistosoma spp* have not been tested in this product, and their associated interference in **parascreen®** is not known.
- 7. Due to limited evidence, the manufacturer does not claim the limit of detection of **parascreen®** for *P.knowlesi* though it is found to detect these Plasmodium species as low as 226 -300 parasites/µl.
- 8. In case of mixed infection (*P.falciparum*, with other malarial species), both, 'Pf' and 'Pan' malaria bands will be positive. Hence, differentiation of infection due to *P.vivax*, *P.ovale* or *P.malariae* cannot be done.
- 9. In *P.falciparum* malaria infection, HRP-2 is not secreted in gametogony stage. Hence, in "Carriers", the HRP-2 band may be absent.










WARRANTY

This product is designed to perform as described on the label and Instructions for use. The manufacturer disclaims any implied warranty of use and sale for any other purpose. In the event of performance changes or product malfunction, please contact manufacturer.

BIBLIOGRAPHY

- 1. Howard, R.J., et al., 1986: Secretion of a Malarial Histidine-rich Protein (Pf. HRP II) from *Plasmodium falciparum*-infected Erythrocytes. J. Cell Biol., 103, 1269-1277.
- 2. Parra, M.E., et al., 1991: Identification of *Plasmodium falciparum* Histidine-Rich Protein 2 in the Plasma of Humans with Malaria. J. Clin. Microbiol., 29, 1629-1634.
- 3. Rodriguez-Del Valle, M., et al., 1991: Detection of Antigens and Antibodies in the Urine of Humans with *Plasmodium falciparum* Malaria. J. Clin. Microbiol., 29, 1236-1242.
- 4. Piper, R. C., et al., (1999) Immuno-capture diagnostic assays for malaria utilizing *Plasmodium* Lactate Dehydrogenase (pLDH) Am. J. Trop. Med. Hyg. 60(1) 109-118.
- 5. Hunte-Cooke A., et al., (1999) Comparison of a Parasite Lactate Dehydrogenase-based Immunochromatographic Antigen Detection assay (OptiMAL®) with Microscopy for the Detection of Malaria Parasites in Human Blood Samples. Am J. Trop Med 60(2). 173-176.
- 6. Quintana M., et al., (1998) Malaria diagnosis by dipstick assay in a Honduran Population with coendemic *Plasmodium falciparum* and *Plasmodium vivax*. Am. J. Trop. Med. Hyg. 59(6) 868-871.
- 7. Palmer, C. J., (1998) Evaluation of OptiMal test for rapid diagnosis of *Plasmodium vivax* and *Plasmodium falciparum*. J. Clin Microbiol. 36(1) 203-206.
- 8. Moody A., et al., (2000) Performance of the OptiMAL™ malaria antigen capture dipstick for malaria diagnosis and treatment monitoring. British Journal of Hematology, 109, 1-5.
- 9. Data on file: Zephyr Biomedicals.

SYMBOL KEYS

|   |  |  |  |
|---|--|--|--|
|  Temperature Limitation            |  Manufacturer                 | <div>DEVICE</div> Device   | <div>EC REP</div> Authorised Representative in the European Community  |
|  Use by                            |  Consult Instructions for use | <div>PIPETTE</div> Disposable Plastic Specimen Applicator  |  |
|  Date of Manufacture               | <div>REF</div> Catalogue Number  | <div>BUF</div> Clearing Buffer   | <div>Xn</div> <div></div> Harmful if swallowed. Do not breathe vapour. If swallowed, seek medical advice immediately and show this container or label. Avoid release to the environment. Refer to special instructions. |
| <div>LOT</div> Batch Number / Lot Number  | <div>IVD</div> In vitro Diagnostic Medical Device  | <div>↑↑</div> This side up   |  |
|  Contains sufficient for <n> tests |  Do not reuse                 |  Do not use if package is damaged |  |



Manufactured by:

**Zephyr Biomedicals**

A Division of Tulip Diagnostics (P) Ltd.

M 46-47, Phase III B, Verna Industrial Estate, Verna, Goa - 403 722, INDIA.

**Regd. Office:** Gitanjali, Tulip Block, Dr. Antonio Do Rego Bagh, Alto Santacruz,

Bambolim Complex P.O., Goa - 403 202, INDIA.

Email address: sales@tulipgroup.com

Tel. : (0832) 2458546, (0832) 2458547

EC REP

CMC Medical Devices & Drugs S.L., Spain.