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**3. Desiccat pouch.**  
4. Package insert.

**STORAGE AND STABILITY**  
The sealed pouches in the test kit may be stored between 4-30°C till the duration of the shelf life as indicated on the pouch / carton. DO NOT FREEZE.

**NOTE**

- For in vitro diagnostic use only. NOT FOR MEDICINAL USE.
- Do not use beyond expiry date.

**SPECIMEN COLLECTION AND PREPARATION**  
Though random urine specimens can be used, first morning urine specimen is preferable as it contains the highest concentration of hCG. Specimens should be collected in clean glass or plastic containers. If testing is not immediate, the urine specimens may be stored at 2-8°C for up to 72 hours. Turbid specimens should be centrifuged or allowed to settle and only the clear supernatant should be used for testing.

**TESTING PROCEDURE AND INTERPRETATION OF RESULTS**

- Bring the sealed pouches to room temperature before testing. Open the pouch and remove the device, dropper and desiccant. Check the colour of the desiccant. It should be blue. If it has turned colourless or pink discard the device and use another device. **Once opened, the device must be used immediately.**
- Refrigerated specimens must be brought to room temperature prior to testing.
- Dispense two drops of urine specimen into the sample well 'S' using the dropper provided.
- At the end of five minutes read the results as follows:

	<b>NEGATIVE:</b> Only one pink/purple colored band appears on the device.
	<b>POSITIVE:</b> Two distinct pink/purple colored bands appear on the device.
	<b>INVALID:</b> The test should be considered invalid if the control band 'C' does not appear. The test is also invalid if only the test band 'T' and no control band 'C' appears. <b>Repeat the test with a new device.</b>

**RESULTS**

- 1) LH calibrator (value : 250 mIU/ml) and dilutions and negative control

LH	250 mIU/ml	125 mIU/ml	60 mIU/ml	30 mIU/ml	15 mIU/ml	7.5 mIU/ml	Contro Neg (n=2)
Results at 5 mins	Pos (55 secs)	Pos (1.06 mins)	Pos (1.16 mins)	Pos (fraca mente)	Negativo	Negativo	

- 2) LH calibrator (value : 550 mIU/ml) and dilution

LH	550 mIU/ml	275 mIU/ml
Results at 5 minutes	Negativo	Negativo

**REFERENCES**

- (1) Batzer, F. R., Hormonal evaluation of early pregnancy, Fertility and Sterility, July 1980), 34, 1. (2) Thompson, R. J., Jackson, A. P., Langlois, N., 1986, Circulating antibodies to mouse monoclonal immunoglobulins in Normal subjects- incidence, species, specificity and effects on a two-site assay for creatine kinase-MB isoenzyme, Clin. Chem. 32, 476-481. (3) Data on file : Viola Diagnostic Systems.

**INTRODUCTION**  
**Tina™ hCG** one step pregnancy test is a rapid, qualitative, two site sandwich immunoassay for the determination of human chorionic gonadotropin (hCG), a marker for pregnancy, in urine specimens.

**SUMMARY**  
Human chorionic gonadotropin (hCG), a glycoprotein hormone secreted by viable placental tissue during pregnancy, is excreted in urine approximately 20 days after the last menstrual period. The levels of hCG rise rapidly reaching peak levels after 60-80 days. The appearance of hCG in urine soon after conception and its rapid rise in concentration makes it an ideal marker for the early detection and confirmation of pregnancy. However elevated hCG levels are frequently associated with trophoblastic and non-trophoblastic neoplasms and hence these conditions should be considered before a diagnosis of pregnancy can be made.

**Tina™ hCG** one step pregnancy test detects the presence of hCG in urine specimens, qualitatively, at concentrations as low as ~25 mIU/ml. Concentration of about 100 mIU/ml of hCG are reached by the first day of the missed menstrual period in normal pregnancy. Thus **Tina™ hCG** one step pregnancy test is able to detect pregnancy at very early stage.

**LIMITATIONS OF THE TEST**

- A number of conditions other than pregnancy including trophoblastic and non-trophoblastic neoplasms such as hydatidiform mole, choriocarcinoma etc., cause elevated levels of hCG. Such clinical conditions must be ruled out before a diagnosis of pregnancy can be made.
- Highly dilute urine specimens and specimens from very early pregnancy may not contain representative levels of hCG. If pregnancy is still suspected, repeat the test with first morning urine after 48-72 hours.
- As with all diagnostic tests, the results must be correlated with clinical findings.

**PRINCIPLE**  
**Tina™ hCG** one step pregnancy test utilizes the principle of IMMUNOCROMATOGRAFIA, a unique two site immunoassay on a membrane. As the test sample flows through the membrane assembly of the device, the colored anti-hCG-coloidal gold conjugate complexes with the hCG in the sample. This complex moves further on the membrane to the test region where it is immobilized by the anti-hCG coated on the membrane leading to formation of a pink/purple colored band which confirms a positive result. Absence of this colored band in the test region indicates a negative test result. The unreacted conjugate and the unbound complex if any move further on the membrane and are subsequently immobilized by the anti-rabbit antibodies coated on the membrane at the control region, forming a pink / purple band. This control band serves to validate the test results.

**PERFORMANCE CHARACTERISTICS**

- Sensitivity:**  
**Tina™ hCG** one step pregnancy test detects the presence of hCG in urine specimens, qualitatively, at concentrations as low as ~25 mIU/ml. Concentration of about 100 mIU/ml of hCG are reached by the first day of the missed menstrual period in normal pregnancy. Thus **Tina™ hCG** one step pregnancy test is able to detect pregnancy at very early stage.
- Specificity:**  
Healthy men and healthy non-pregnant women do not have detectable levels of hCG by **Tina™ hCG** one step pregnancy test.
- Accuracy:**  
The results obtained by **Tina™ hCG** one step pregnancy test correlated very well when run in parallel with other commercially available tests for pregnancy, using known positive and negative specimens.

**CONSULT INSTRUCTIONS FOR USE**

**Manufacturer**

**This side up**

**REF**

**LOT**

**Do not reuse**

**Read results at the end of five minutes**

**IVD**

**Manufactured by:**  
**Viola Diagnostic Systems**  
A Division of Tulip Diagnostics (P) Ltd.  
Plot No. 33, Sector-3, I.I.E. SIDCUL, Pantharganj, U.S. Nagar, Uttarakhand - 263 153, INDIA.  
Regd. Office: Gianjali, Tulla Block, Dr. Antonio Do Rego Bagh, Alto Santacruz, Bambolim Complex P.O., Goa - 403 202, INDIA.

**EC REP**

CMC Medical Devices & Drugs S.L., Spain.