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Tina™ hCG ONE STEP PREGNANCY TEST (Dipstick)

INTRODUÇÃO

Tina™ hCG one step pregnancy test is a rapid, qualitative, two site sandwich immunoassay for the determination of human chorionic gonadotropin (hCG), a marker for pregnancy, in urine specimens.

SUMMARY

Human chorionic gonadotropin (hCG), a glycoprotein hormone secreted by viable placental tissue during pregnancy, is excreted in urine approximately 20 days after the last menstrual period. The levels of hCG rise rapidly reaching peak levels after 60-80 days. The appearance of hCG in urine soon after conception and its rapid rise in concentration makes it an ideal marker for the early detection and confirmation of pregnancy. However elevated hCG levels are frequently associated with trophoblastic and non-trophoblastic neoplasms and hence these conditions should be considered before a diagnosis of pregnancy can be made.

Tina™ hCG one step pregnancy test detects the presence of hCG in urine specimens, qualitatively, at concentrations as low as ~25 mIU/ml in less than five minutes.

PRINCIPLE

Tina™ hCG one step pregnancy test utilizes the principle of IMMUNOCHEMOTROPHOGRAPHY, a unique two site immunoassay on a membrane. As the test sample flows through the membrane assembly of the dipstick, the colored anti-hCG-colloidal gold conjugate complexes with the hCG in the sample. This complex moves further on the membrane to the test region where it is immobilized by the anti-hCG coated on the membrane leading to formation of a pink/purple colored band which confirms a positive test result. Absence of this colored band in the test region indicates a negative test result. The unreacted conjugate and the unbound complex if any move further on the membrane and are subsequently immobilized by the anti-rabbit antibodies coated on the membrane at the control region, forming a pink / purple band. This control band serves to validate the test results.

REAGENTS AND MATERIALS SUPPLIED

Each individual foil packed pouch contains:

1. Dipstick: Membrane assembly predispensed with anti-hCG

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NOTA

1. For in vitro diagnostic use only. NOT FOR MEDICAL USE.
2. Do not use beyond expiry date.

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Tina™ hCG ONE STEP PREGNANCY TEST (Dipstick)

INTRODUCTION

Tina™ hCG one step pregnancy test is a rapid, qualitative, two site sandwich immunoassay for the determination of human chorionic gonadotropin (hCG), a marker for pregnancy, in urine specimens.

SUMMARY

Human chorionic gonadotropin (hCG), a glycoprotein hormone secreted by viable placental tissue during pregnancy, is excreted in urine approximately 20 days after the last menstrual period. The levels of hCG rise rapidly reaching peak levels after 60-80 days. The appearance of hCG in urine soon after conception and its rapid rise in concentration makes it an ideal marker for the early detection and confirmation of pregnancy. However elevated hCG levels are frequently associated with trophoblastic and non-trophoblastic neoplasms and hence these conditions should be considered before a diagnosis of pregnancy can be made.

Tina™ hCG one step pregnancy test detects the presence of hCG in urine specimens, qualitatively, at concentrations as low as ~25 mIU/ml in less than five minutes.

PRINCIPLE

Tina™ hCG one step pregnancy test utilizes the principle of IMMUNOCHEMOTROPHOGRAPHY, a unique two site immunoassay on a membrane. As the test sample flows through the membrane assembly of the dipstick, the colored anti-hCG-colloidal gold conjugate complexes with the hCG in the sample. This complex moves further on the membrane to the test region where it is immobilized by the anti-hCG coated on the membrane leading to formation of a pink/purple colored band which confirms a positive test result. Absence of this colored band in the test region indicates a negative test result. The unreacted conjugate and the unbound complex if any move further on the membrane and are subsequently immobilized by the anti-rabbit antibodies coated on the membrane at the control region, forming a pink / purple band. This control band serves to validate the test results.

ACCURACY:

The results obtained by Tina™ hCG one step pregnancy test correlated very well when run in parallel with other commercially available tests for pregnancy, using known positive and negative specimens.

antiserum-colloidal gold conjugate, rabbit IgG colloidal gold conjugate and anti-hCG antisera and anti-rabbit antisera at the respective regions.

STORAGE AND STABILITY

1. Store the product at 4-30°C, away from heat and moisture, preferably in a cool dark place. The foil packed dipsticks in unopened bottles are stable up to the expiry date as mentioned on the product label.
2. Once the bottle seal is opened, the unopened foil packed dipsticks must be used immediately as directed.
3. The desiccants provided in the foil pack should be reduced moisture in the bottle while packing, storage and transit.
4. Once the bottle seal is opened, the desiccants are dysfunctional, however the foil packed strips are stable for a 60 day period.

NOTE

- 1