



<b>ADAZYME™-LS</b>	<b>Enzymatic Assay For ADA Determination</b>												
<b>Application</b>	For the Determination of Adenosine Deaminase Activity in Serum, Plasma and Biological Fluids by Enzymatic Method												
<b>Principle</b>	The ADA assay is based on the enzymatic deamination of adenosine to inosine which is converted to hypoxanthine by purine nucleoside phosphorylase (PNP). Hypoxanthine is then converted to uric acid and hydrogen peroxide (H <sub>2</sub> O <sub>2</sub> ) by xanthine oxidase (XOD). H <sub>2</sub> O <sub>2</sub> is further reacted with N-Ethyl-N-(2-hydroxy-3-sulfopropyl)-3-methylaniline (EHSPT) and 4-amino antipyrine (4-AA) in the presence of peroxidase (POD) to generate a Quinone dye which is monitored in a kinetic manner												
<b>Specificity</b>	<table border="0" style="width: 100%;"> <tr> <td style="width: 40%;">Adenosine + H<sub>2</sub>O</td> <td style="width: 10%; text-align: center;">— ADA —&gt;</td> <td style="width: 50%;">Inosine + NH<sub>3</sub></td> </tr> <tr> <td>Inosine + Pi</td> <td style="text-align: center;">— PNP —&gt;</td> <td>Hypoxanthine + Ribose-1-phosphate</td> </tr> <tr> <td>Hypoxanthine + H<sub>2</sub>O + O<sub>2</sub></td> <td style="text-align: center;">— XOD —&gt;</td> <td>Uric acid + H<sub>2</sub>O<sub>2</sub></td> </tr> <tr> <td>H<sub>2</sub>O<sub>2</sub> + 4-AA + EHSPT</td> <td style="text-align: center;">— POD —&gt;</td> <td>H<sub>2</sub>O + Quinone dye (λ<sub>max</sub> 546nm)</td> </tr> </table> <p>One unit of ADA is defined as the amount of ADA that generates one μmole of inosine from adenosine per min at 37°C.</p>	Adenosine + H <sub>2</sub> O	— ADA —>	Inosine + NH <sub>3</sub>	Inosine + Pi	— PNP —>	Hypoxanthine + Ribose-1-phosphate	Hypoxanthine + H <sub>2</sub> O + O <sub>2</sub>	— XOD —>	Uric acid + H <sub>2</sub> O <sub>2</sub>	H <sub>2</sub> O <sub>2</sub> + 4-AA + EHSPT	— POD —>	H <sub>2</sub> O + Quinone dye (λ <sub>max</sub> 546nm)
Adenosine + H <sub>2</sub> O	— ADA —>	Inosine + NH <sub>3</sub>											
Inosine + Pi	— PNP —>	Hypoxanthine + Ribose-1-phosphate											
Hypoxanthine + H <sub>2</sub> O + O <sub>2</sub>	— XOD —>	Uric acid + H <sub>2</sub> O <sub>2</sub>											
H <sub>2</sub> O <sub>2</sub> + 4-AA + EHSPT	— POD —>	H <sub>2</sub> O + Quinone dye (λ <sub>max</sub> 546nm)											
<b>Reagent System</b>	a) R1 - ADAZYME -LS: Enzyme Reagent, ready to use. b) R2 - ADAZYME -LS : Starter Reagent, ready to use c) C - ADAZYME -LS : Calibrator (Lyophilized)												
<b>Test duration</b>	11 mins												
<b>Storage / Stability</b>	2-8°C, 18 months												
<b>In Use Stability</b>	Reconstituted Calibrator is stable for 15 days at 2-8°C.												

Presentation	Pack	Cat.No.
Adazyme™ LS	10 ml	1102310010
	25 ml	1102310025