

## Hemospot® Kit Standard Guaiac Method

### Intended Use:

Hemospot® test is useful in the detection of bleeding caused by gastrointestinal disorders such as colitis, polyps, diverticulitis and hookworm infestation. Fecal occult blood tests are recommended for use as an aid to routine physical examinations. Routine hospital testing, screening for gastrointestinal bleeding from any source including colorectal cancer. Hemospot® kit uses standard guaiac method to determine occult blood in stool.

### Hemospot® Kit components:

<b>Test Cards</b>	Hemospot® test cards consisting of a filter paper impregnated with the guaiac resin (reactive surface)
<b>Developer solution</b>	Dropper bottle containing developer solution, ready to use
<b>Positive control</b>	Ready to use
<b>Other Accessories</b>	Sample applicator and package insert

### Linearity:

The reaction area turns blue if occult blood is present in the sample

FEATURES	BENEFITS
<b>Convenient</b>	No sample preparation required
<b>Easy procedure and result interpretation</b>	2-step rapid; add stool sample, dry & add developer solution No instrumentation required

Storage / Stability	Temperature	Duration
<b>Unopened kit</b>	15-25°C	24 Months
<b>Opened kit</b>	15-25°C	24 Months

Available Pack Sizes	
25 Tests	50 Tests