

Quantia® Cystatin C Control Set For QC

Intended Use: _____

For monitoring the precision of laboratory testing procedure for Cystatin C quantitation, for validating the calibration and performance of Quantia® Cystatin C reagents.

Quantia® Cystatin C Control Set:

Quantia® Cystatin C control (Level I, Level II)	2 level ready to use controls
Other Accessories	Package Insert

FEATURES	BENEFITS
Liquid Stable Controls	Convenient Eliminates possibilities of reconstitution error Better recovery of values
A set of two level controls	Ensure complete QC at higher and lower level of Cystatin C assays

Storage/ Stability	Temperature	Duration
Unopened	2-8°C	12 months
Open vial stability	2-8°C	7 days when stored tightly capped

Available Pack Sizes
2 x 0.5 ml