

Quantia[®] HbA1c Calibrator Set

For calibration of HbA1c assays in human blood

Intended Use:

Quantia[®] HbA1c calibrator set is intended for use as a calibrator (standard) set for preparation of standard curves for quantitative determination of % HbA1c.

Quantia[®] HbA1c calibrator set components:

Quantia[®] HbA1c calibrator (Level I, Level II, Level III)	Lyophilized calibrators with different HbA1c concentrations
Other Accessories	Package Insert

FEATURES	BENEFITS
Calibrator traceable to DCCT	Results aligned to NGSP/IFCC Good inter-method agreement in results
Set of three level calibrator	Accurate values
Freeze dried reagent	Aids in minimizing microbial contamination Improves reagent stability

Storage/ Stability	Temperature	Duration
Unopened	2-8°C	12 months
Reconstituted reagent	2-8°C	7 days when stored tightly capped

Available Pack Sizes
3 x 0.2 ml