

Redgen® Toluidine Red Unheated Serum Test (TRUST) for rapid serodiagnosis of syphilis

Intended Use: _____

Redgen®, a Toluidine Red Unheated Serum Test (TRUST) card test for the detection of Anti-Lipoidal antibodies in serum and plasma during syphilis.

Redgen® Kit components:

Reagent	Microfine toluidine red dye suspension coated with lipid complexes
Positive Control	Reactive with the reagent
Negative Control	Non Reactive with the reagent
Slide	Disposable slide with eight reaction circles
Other Accessories	Dropper, pipette, mixing sticks, rubber teats and Pack Insert

Performance Characteristics:

Specificity: Lipid complex reactive with Anti-Lipoidal antibodies

FEATURES	BENEFITS
Specially designed anti - complement agent removes complement interference	No false reactions due to complements
Coated disposable cards provide optimum reaction environment	Unmatched readability
Reagents standardized against WHO IRP	Reliable performance of kit
Special pack design	Saves refrigerator space

Storage / Stability	Temperature	Duration
Unopened kit	2-8°C	15 Months
Opened kit	2-8°C	15 Months

Available Pack Sizes
75 Tests