

## Multidet Multipurpose neutral laboratory detergent

**Intended Use:** \_\_\_\_\_

Multipurpose neutral laboratory detergent for cleaning sensitive equipment / parts, glass precision equipments, glassware.

**Multidet:**

<b>Multidet</b>	Concentrated neutral detergent
-----------------	--------------------------------

FEATURES	BENEFITS
Appropriate	Cleans the glassware and other laboratory equipments with minimum residue
Flexible usage	Can be used with different glassware, precision equipments, etc.
Convenient pack size	Suitable for all users

<b>Suggested concentration</b>	2.5%
--------------------------------	------

Storage / Stability	Temperature	Duration
Unopened kit	25-30°C	24 Months
Opened kit	25-30°C	24 Months

Available Pack Sizes
1000 ml