

**Gravichcek[®]
Rapid one step test dipstick for hCG**

Intended Use: _____

Gravichcek[®], a rapid one step test dipstick for the detection of hCG in urine/serum during pregnancy.

Gravichcek[®] Dipstick kit components:

Dipstick	Membrane assembly pre-dispensed with monoclonal anti-hCG antisera - colloidal gold conjugate, rabbit IgG colloidal gold conjugate monoclonal anti-hCG antiserum at the 'Test region' and anti-rabbit anti serum at the 'Control region'
Pipette	Disposable plastic sample dropper
Other Accessories	Desiccant pouch and Package Insert

Performance Characteristics:

Sensitivity: ~ 10 mIU/ml of hCG	Specificity: 100% specific to β hCG
--	--

FEATURES	BENEFITS
Sensitivity adjusted to detect > 10 mIU/ml of hCG	Optimum sensitivity for safe and confident reporting of results Facilitates early pregnancy detection
Specific for beta sub unit of hCG molecule	No cross-reactivity with other hormones
Dual sample options: Urine / Serum	Versatile / Offers specimen flexibility Routine pregnancy testing using urine specimens Enables diagnosis of problem pregnancies
Additional benefits of serum testing	Accurate, reliable and confident reporting of results even in problem pregnancies and bad obstetrics cases

Storage / Stability	Temperature	Duration
Unopened kit	4-30°C	24 Months
Opened kit	4-30°C	< 1 Hour

Available Pack Sizes
25 Tests