

HbChek

Smart Hemoglobin Monitoring System



KEY FEATURES

- Gold Standard in Point of Care Testing
- State of The Art Microfluid Technology
- High Accuracy – CV < 3% As Compared with Gold Standard Method
- Simple to Use Requires Very Low Maintenance
- Test Sample – Only 8µl Blood Sample
- Smart Connectivity – Connect with Mobile App Through Smart Bluetooth
- Easy to View Result – Large LCD Display
- 1000 Result Storage Memory.

Technical Specifications

Test Principle	Transmittance/ Absorbance Spectrophotometry
Measuring Range	0 to 25 g/dL
Sample Type	Capillary Blood or Venous Blood
Sample Volume	~ 8µl
Precision	CV < 3 %
Reading Time	25 Seconds
Power	3 X AAA Alkaline Batteries
Weight	~ 250gms.
Operating Temperature	0 – 50 ° C
Connectivity	Bluetooth Smart BLE
Calibration	Factory Calibrated, No Need to Calibrate in Field

LIST OF COMPONENTS SUPPLIED WITH THE INSTRUMENT

Components	Quantity
Strivettes	20 No.
Lancing Device	1 No.
Lancet	20 No.
Alcohol Swabs	20 No.
Tray Cleaning Tool	1 No.
Manual	1 No.
Batteries	3 No.

ACCESSORIES/CONSUMABLES FOR USE WITH VEINSPY

Accessories/Consumables	Cat No.	Pack Size
Disposable Strivettes	847HM00002	20 Test