

Microxpress®



Pouch Warmer PW-004

KEY FEATURES

- No requirement of water bath.
- Proper melting of media ensures correct results in microbiology.
- Calibrated & standardized for Instaprep media
- Seven segment LED display with front panel keys.
- Timer with alarm set to 15 minutes.

TECHNICAL SPECIFICATIONS

Principle	Incubation
Pouch capacity	4
Warm up time	Maximum 30 minutes
Working temperature	90° C
Operating conditions	Temperature- Ambient 25° - 30° C
Power supply	230 V / 50 Hz
Dimensions	25.5 cms x 19.5 cms x 18 cms
Net Weight	5.3 Kgs
Cat No.	T10025

LIST OF STANDARD ACCESSORIES SUPPLIED ALONG WITH THE INSTRUMENT

Components	Qty
Fuses (5A)	2 Nos.
User Manual	1 No.