



matrixTM
—GEL•SYSTEM

Card Warmer

CW-2400

KEY FEATURES

- Seven segment LED display with front panel keys.
- Working temperature 37°C (set temp.)
- Warm up time maximum 30 minutes.
- Timer with alarm (15 minutes).

TECHNICAL SPECIFICATIONS

Principle	Incubation
Card capacity	24 cards
Operating conditions	Temperature- Ambient 25° - 40°C
Power Supply	230V, 50 Hz.
Dimensions	29 cms x 25.5 cms x 14.5 cms
Net Weight	5.1 Kgs
Cat. Nos.	821CW002400

LIST OF STANDARD ACCESSORIES SUPPLIED ALONG WITH THE INSTRUMENT

Components	Qty
Power Cord	1 No.
User Manual	1 No.