



**matrix**<sup>TM</sup>  
—GEL•SYSTEM

## Card Warmer

CW-2400

### KEY FEATURES

- Seven segment LED display with front panel keys.
- Working temperature 37°C (set temp.)
- Warm up time maximum 30 minutes.
- Timer with alarm (15 minutes).

### TECHNICAL SPECIFICATIONS

<b>Principle</b>	Incubation
<b>Card capacity</b>	24 cards
<b>Operating conditions</b>	Temperature- Ambient 25° - 40°C
<b>Power Supply</b>	230V, 50 Hz.
<b>Dimensions</b>	29 cms x 25.5 cms x 14.5 cms
<b>Net Weight</b>	5.1 Kgs
<b>Cat. Nos.</b>	821CW002400

### LIST OF STANDARD ACCESSORIES SUPPLIED ALONG WITH THE INSTRUMENT

Components	Qty
Power Cord	1 No.
User Manual	1 No.