



PYROGUN™

NON-CONTACT THERMOMETER

KEY FEATURES

- It can measure human body temperature & is quite helpful for monitoring temperature of pediatric patient as well as agitated patients.
- Conveniently measures temperature of patients in hospital OPDs, as it measure temperature accurately in less than one second.
- Liquid temperature such as, milk for infants, bath water can be measured accurately.
- Surface temperature of floor, incubators, food can be measured whenever required.

TECHNICAL SPECIFICATIONS

Measuring Range	Body Mode : 90°F-109°F(32°C - 42°C) Surface Mode : 32°F-140°F(10°C - 60°C)
Accuracy	±0.1°C (0.3°F)
Measuring Distance	5 -15 Cm (2" - 6")
Measuring time	Less than 1 Second
3 Color backlight (Color Alarm) <i>Applicable for "Body mode"</i>	Green back Light(Normal): Upto 99.1°F (37.2°C) Orange back light (Slight Fever) From 99.2°F to 100.4°F (37.2°C - 38.°C) Red back light (High Fever) Above 100.4°F (38°C)
Data Storage	Can Store 34 recorded temperature
Auto power off	7 seconds
Battery	3 V ; AA ; Alkaline Battery
Weight	122 gm
Size	145 mm X 80mm X 40 mm