



## The Hand-Held Vein Finding Device

### KEY FEATURES

- Reduces multiple needle pricks & improves site selection for needle insertion.
- No treatment delay due to difficult venous access
- Reduces patient pain due to needle insertion in muscle/ tingling sensation due to hitting of needle in nerve.
- Increases patient comfort & satisfaction. Avoids hematoma.

### TECHNICAL SPECIFICATIONS

<b>Number of LED's</b>	22 Orange, 8 Red
<b>Power ON</b>	Indication by Green LED
<b>Low power</b>	Indication by Red LED
<b>Functional opening</b>	17mm
<b>Viewing diameter</b>	33mm
<b>Rechargeable Battery</b>	Lithium Ion, 3.7V, 1150mAh
<b>Battery charger</b>	110-240 volts AC input, 5V DC output 750mA
<b>Battery Backup</b>	1½ hrs. Continuous usage.
<b>Dimensions</b>	L 125 mm; B 65 mm; H 22.7 mm
<b>Net Weight</b>	90 gms.

### CONSUMABLES AVAILABLE ON ORDER

<b>Components</b>	<b>Cat No.</b>	<b>Pack Size</b>
Disposable Plastic Sleeves	825DS5000050	50 Nos