

**Dengue NS-1****Rapid test for the detection of NS-1 antigen to Dengue virus in human serum/plasma**

Intended Use: _____

Dengue NS-1 is intended for the detection of NS-1 antigen to Dengue virus in human serum/plasma.

Dengue NS-1 kit components:

Device	Membrane assembly for Dengue NS-1 antigen is predispensed with: <ul style="list-style-type: none"> ● Monoclonal 'anti-dengue NS-1' colloidal gold conjugate ● Monoclonal 'anti-dengue NS-1' antibody at the test window 'T' ● Anti- mouse globulin predispensed at the control window 'C'
Dropper	Disposable plastic dropper
Other Accessories	Desiccant pouch and Package Insert

Performance Characteristics:

Sensitivity: 100%	Specificity: 100%
--------------------------	--------------------------

FEATURES	BENEFITS
Test system designed to detect dengue NS-1 antigen	Detects dengue infection on the first day of appearance of clinical symptoms.
No sample running buffer required	Easy, user friendly procedure
Sensitivity 100% and specificity is 100%	Good sensitivity and specificity
Pack size 10 test with a shelf life of 24 months	Convenient for small laboratories also
Storage at 4-30°C	Well suited for field conditions, on-site diagnosis, flexibility of remote distribution.

Storage / Stability	Temperature	Duration
Unopened kit	4-30°C	24 Months
Opened pouch	Room temperature	Use immediately

Available Pack Sizes
10 Tests