



SEROCHECK Tp™
Rapid test for Syphilis in human serum / plasma.

Intended Use: _____

Serocheck Tp™ – Dipstick, a rapid, qualitative, two site double antigen sandwich Immunoassay (modified TPHA) for the detection of *Treponema* specific IgG and IgM antibodies in human serum or plasma specimens.

Serocheck Tp™ – Dipstick Kit Components:

Dipstick	Each individual pouch contain Membrane assembly pre-dispensed with: <ul style="list-style-type: none"> ● Recombinant <i>Treponema pallidum</i> antigen colloidal gold conjugate. ● Recombinant <i>Treponema pallidum</i> antigen at 'Test' region. ● Agglutinating sera for rabbit globulin coated at 'Control' region.
Other accessories	Desiccant pouch, package insert and dipstick holder.

Performance Characteristics:

Sensitivity: 100%	Specificity: 98%
--------------------------	-------------------------

FEATURES	BENEFITS
Highly sensitive <i>Treponema</i> specific recombinant antigens used.	<ul style="list-style-type: none"> ● Early detection. ● Ideal for screening and confirmation
Double Antigen Sandwich system.	<ul style="list-style-type: none"> ● Detects IgM & IgG class of antibodies. ● Detects Syphilis infection at every stage of the disease.
Evaluated with Syphilis ACCUSET Performance Panel (0820-0214) and Syphilis SEROCONVERSION Panel (0615-0017/PSS901-1.2) obtained from Sera care Life Sciences	<ul style="list-style-type: none"> ● Trustworthy Performance. ● Reliability.
No specimen Dilution required	<ul style="list-style-type: none"> ● No variability ● Less errors. ● No loss in sensitivity
Unique design - Safety Dipstick Holder	<ul style="list-style-type: none"> ● Prevents contamination of the testing surface. ● High safety for the user. ● Easy to discard.

Storage / Stability	Temperature	Duration
Unopened kit / pouch	4-30°C	24 Months
Opened pouch	4-30°C	Use immediately

Available Pack Sizes
25 Tests