



Zephyr Biomedicals

**Leishcheck™**  
**Rapid test for detection of Visceral Leishmaniasis antibody in human serum/plasma**

**Intended Use:** \_\_\_\_\_

Leishcheck™ is intended for the presumptive diagnosis of Visceral Leishmaniasis (VL).

**Leishcheck™ kit components:**

<b>Device</b>	Membrane assembly predisposed with Protein A- colloidal gold conjugate, recombinant K39 antigen (VL specific) and reagent control coated at the respective regions
<b>Buffer</b>	Ready to use buffer in a dropper bottle
<b>Other Accessories</b>	Desiccant pouch and Package Insert

**Performance Characteristics:**

<b>Sensitivity: 100%</b>	<b>Specificity: 99.5%</b>
--------------------------	---------------------------

<b>FEATURES</b>	<b>BENEFITS</b>
<b>Recombinant K39(rK39) antigen used</b>	Detects specific antibodies to Leishmaniasis
<b>Convenient for all class of laboratories</b>	Pack size 10 tests with a shelf life of 24 months
<b>Storage at 4-30°C</b>	Well suited for field conditions, on-site diagnosis, flexibility of remote distribution

<b>Storage / Stability</b>	<b>Temperature</b>	<b>Duration</b>
<b>Unopened kit</b>	4-30°C	24 Months
<b>Opened device</b>	R.T.	Use immediately

<b>Available Pack Sizes</b>
10 Tests