

CE

TinaTM hCG

One Step Pregnancy Test



GLOBALLY TRUSTED

TinaTM hCG


Simple • Rapid • Reliable

How to do the test

Step 1 Collect urine



Step 2 Add 2 drops to test cassette



Step 3 Wait for 5 mins and read results



FEATURES	BENEFITS
Sensitivity calibrated to detect ~25mIU/ml of hCG	Ideal for early pregnancy detection Exclude normal non pregnant individuals Facilitates safe and confident reporting of results
Simple and Accurate	One step test without requirement of skill Correlates very well with other commercially reliable kits
Highly specific	100% specificity for β hCG No cross reactivity with LH and FSH.
Built in control	Automated validation of the test performance
Composite Pack	Individually pouched cassettes ensuring moisture-free protection.
CE certified	Extremely Reliable performance
Results within 5 minutes	Quick and rapid
Storage at 4-30°C	Facilitates remote distribution and on the spot testing



For further information contact:

T TULIP DIAGNOSTICS (P) LTD

Gitanjali, Tulip Block, Dr. Antonio Do Rego Bagh, Alto Santacruz, Bambolim Complex Post Office, Goa - 403202, INDIA.
Tel.: +91 832 2458546-50. E-mail: sales@tulipgroup.com Website: www.easyhometests.in

orchid



Microexpress

Coral

Clinical Systems

BioShields

Viola

There is no alternative to a better pregnancy test