



### Agtrol® Starter pack for preparing Coombs control cells

Intended Use: \_\_\_\_\_

Agtrol® is an appropriately diluted monoclonal Anti-D IgG and red cell preserving solution starter pack for preparation of Coombs control cells in situ. The usage of Coombs cells is advocated for functional validation of anti-human globulin reagent and procedural validation of tests employing anti-human globulin reagent.

Agtrol® Kit components:

Anti-D (IgG) monoclonal antibody reagent	Ready to use reagent standardized pre-diluted solution of Anti-D (IgG) monoclonal antibody reagent
R.B.C. preserving solution	Ready to use solution, red blood cell preserving solution for serological applications
Other Accessories	Package Insert

Performance Characteristics:

Standardization: Follows AABB and FDA guidelines

FEATURES	BENEFITS
Ready to use solutions	No manual preparation of reagent and preserving solution required
Dual Purpose	Functional validation of anti-human globulin reagent Procedural validation of tests employing anti-human globulin reagent

Storage/Stability	Temperature	Duration
Unopened kit	2-8°C	24 Months
Opened kit	2-8°C	24 Months

Available Pack Sizes
10 ml