

**Erybank[®] Rh-hr Control****Negative control for Rh₀ (D) typing (slide and tube test)****Intended Use:** _____Erybank[®] Rh-hr control is used to create a negative control parallel to Anti-D testing using Anti-D (Rh₀) IgM and Anti-D (Rh₀) IgM+IgG.**Erybank[®] Rh-hr Control Kit components:**

Erybank[®] Rh-hr Control	Ready to use reagent that contains the same ingredients as Anti-D reagent without the specific antibodies.
Other Accessories	Package Insert

Performance Characteristics:

Total Protein: ~ 2% (g/dl)	Standardization: Follows FDA guidelines
-----------------------------------	--

FEATURES	BENEFITS
Contains all ingredients in appropriate protein concentration as same of Anti-D reagent except the specific antibodies.	Ideal reagent for Negative control in Rh ₀ (D) typing in Direct Anti-D testing
No cross reactivity with RBC's	

Storage / Stability	Temperature	Duration
Unopened kit	2-8°C	24 Months
Opened kit	2-8°C	24 Months

Available Pack Sizes	
5 ml	10 ml