

**Erybank<sup>®</sup> Rh-hr Control****Negative control for Rh<sub>0</sub> (D) typing (slide and tube test)****Intended Use:** \_\_\_\_\_Erybank<sup>®</sup> Rh-hr control is used to create a negative control parallel to Anti-D testing using Anti-D (Rh<sub>0</sub>) IgM and Anti-D (Rh<sub>0</sub>) IgM+IgG.**Erybank<sup>®</sup> Rh-hr Control Kit components:**

<b>Erybank<sup>®</sup> Rh-hr Control</b>	Ready to use reagent that contains the same ingredients as Anti-D reagent without the specific antibodies.
<b>Other Accessories</b>	Package Insert

**Performance Characteristics:**

<b>Total Protein:</b> ~ 2% (g/dl)	<b>Standardization:</b> Follows FDA guidelines
-----------------------------------	--

FEATURES	BENEFITS
Contains all ingredients in appropriate protein concentration as same of Anti-D reagent except the specific antibodies.	Ideal reagent for Negative control in Rh <sub>0</sub> (D) typing in Direct Anti-D testing
No cross reactivity with RBC's	

Storage / Stability	Temperature	Duration
Unopened kit	2-8°C	24 Months
Opened kit	2-8°C	24 Months

Available Pack Sizes
5 ml