



## Erygen-AC

### Ready to use Coombs control cells

**Intended Use:**

Erygen-AC, are ready to use Coombs control cells. It is a panel of 3% reagent red blood cells for quality control of anti-human globulin procedures. These erythrocytes, sensitized with Anti-D (IgG) monoclonal (Coombs control cells) are used as positive controls in anti-human globulin testing.

**Erygen-AC Kit components:**

<b>Erygen-AC</b>	Panel of 3% RBCs of group Rho (D) positive cells sensitized with Anti-D (IgG) monoclonal
<b>Other Accessories</b>	Package Insert

**Performance Characteristics:**

<b>Specificity:</b> Standardized to provide positive agglutination control for anti-human globulins containing IgG
<b>Standardization:</b> Follows AABB and FDA guidelines

FEATURES	BENEFITS
<b>Convenient</b>	Available in a hassle-free, easy to use and ready prepared solution Components in precise quantities ensure maximum benefits during quality control procedures
<b>Single reagent pack</b>	Having a shelf life of 50 days and with added preservatives to maintain cell integrity

Storage / Stability	Temperature	Duration
<b>Unopened kit</b>	2-8°C	50 days
<b>Opened kit</b>	2-8°C	50 days

Available Pack Sizes
2 ml