



Matrix™ ABO Subgrouping Card Anti-A₁

Intended Use: _____

Matrix™ ABO Subgrouping Card Anti-A₁ is used for confirmation of A₁ subgroup of A.

Matrix™ ABO Subgrouping Card Anti-A₁ kit components:

Matrix™ ABO Subgrouping Card Anti-A₁	Matrix™ ABO Subgrouping Card Anti-A ₁ contains Anti-A ₁ lectin
Configuration	6 x Anti-A ₁ / Six profiles per card
Other Accessories	Card fitment and Package Insert

Performance Characteristics:

Standardization: Follows FDA guidelines
--

FEATURES	BENEFITS
Reliable and consistent reactivity only with A₁ and A₁B cells	Optimized to exclude A ₂ reactivity and weaker subgroups Calibrated to pick up only A ₁ and A ₁ B cells
Stable End Reactions	Read results hours later Can be saved for QA or peer review Can be photocopied for documentation
Objective Interpretation of results	Easy Interpretation of weak reactions Relatively skill free
Small sample size	More flexible More convenient
Indigenously Manufactured	Easy availability Long shelf life Another first from Tulip
Tulip's Technical Support	Dependability

Storage / Stability	Temperature	Duration
Unopened kit	4-25°C	12 Months
Opened card	---	Use immediately

Available Pack Sizes
24 Cards