



Matrix™ AHG (Coombs) Test Card

Intended Use: _____

Matrix™ AHG (Coombs) Test Card is used for direct coombs test (DAT) and indirect coombs tests (IAT) including antibody screening, antibody identification and compatibility testing. Contains polyspecific coombs reagent.

Matrix™ AHG (Coombs) Test Card kit components:

Matrix™ AHG (Coombs) Test Card	Matrix™ Anti-Human Globulin card contains six microtubes prefilled with gel in a suitable buffer containing Anti-human IgG and Monoclonal Anti-C ₃ d. The gel contains a suitable preservative
Configuration	6 x AHG / Profiles per card depends on test parameter
Other Accessories	Card fitment and Package Insert

Performance Characteristics:

Standardization: Follows FDA guidelines
--

FEATURES	BENEFITS
Stable End Reactions	Read results hours later Can be saved for QA or peer review Can be photocopied for documentation
Objective Interpretation of results	Easy Interpretation of weak reactions Relatively skill free
Enhanced Sensitivity	Better detectability especially in antibody screen and cross match
Eliminates cell wash step in IAT	Less skill and technique demands Saves time
Small sample size	More flexible More convenient
Indigenously Manufactured	Easy availability Long shelf life Another first from Tulip
Tulip's Technical Support	Dependability

Storage / Stability	Temperature	Duration
Unopened kit	4-25°C	12 Months
Opened card	---	Use immediately

Available Pack Sizes
24 Cards