



Matrix™ Coombs and Neutral Gel Card

Intended Use: _____

Matrix™ Coombs and Neutral Gel Card is a combination of Coombs and neutral gel can be used for multiple applications such as antibody screening in AHG and Enzyme phase or Compatibility testing in AHG and Enzyme/saline phases.

Matrix™ Coombs and Neutral Gel Card kit components:

Matrix™ Coombs and Neutral Gel Card	Matrix™ Coombs and Neutral Gel Card contains polyspecific Coombs (Anti-IgG+ monoclonal C3d) and Neutral gel
Configuration	3 x AHG- 3 x Neutral / Profiles per card depends on test parameter
Other Accessories	Card fitment and Package Insert

Performance Characteristics:

Standardization: Follows FDA guidelines

FEATURES	BENEFITS
Stable End Reactions	Read results hours later Can be saved for QA or peer review Can be photocopied for documentation
Objective Interpretation of results	Easy Interpretation of weak reactions Relatively skill free
Enhanced Sensitivity	Better detectability especially in antibody screen and cross match
Small sample size	More flexible More convenient
Indigenously Manufactured	Easy availability Long shelf life Another first from Tulip
Tulip's Technical Support	Dependability

Storage / Stability	Temperature	Duration
Unopened kit	4-25°C	12 Months
Opened card	---	Use immediately

Available Pack Sizes
24 Cards