



### Matrix™ Coombs Anti-IgG Card

**Intended Use:** \_\_\_\_\_

Matrix™ Coombs Anti-IgG card is used for all coombs test where anti complement activity is not required like direct coombs test and indirect coombs test including antibody screening, antibody identification and compatibility testing. Weak D testing can be done using Agtrol.

**Matrix™ Coombs Anti-IgG Card kit components:**

<b>Matrix™ Coombs Anti-IgG Card</b>	Matrix™ Coombs Anti-IgG card contains six microtubes pre-filled with gel in a suitable buffer containing monospecific Anti-IgG antibodies. The gel contains a suitable preservative
<b>Configuration</b>	6 x Anti-IgG / Profiles per card depends on test parameter
<b>Other Accessories</b>	Card fitment and Package Insert

**Performance Characteristics:**

**Standardization:** Follows FDA guidelines

FEATURES	BENEFITS
<b>Stable End Reactions</b>	Read results hours later Can be saved for QA or peer review Can be photocopied for documentation
<b>Objective Interpretation of results</b>	Easy Interpretation of weak reactions Relatively skill free
<b>Enhanced Sensitivity</b>	Better detectability especially in antibody screen and cross match
<b>Eliminates cell wash step in IAT</b>	Less skill and technique demands Saves time
<b>Small sample size</b>	More flexible More convenient
<b>Indigenously Manufactured</b>	Easy availability Long shelf life Another first from Tulip
<b>Tulip's Technical Support</b>	Dependability

Storage / Stability	Temperature	Duration
Unopened kit	4-25°C	12 Months
Opened card	---	Use immediately

Available Pack Sizes
24 Cards