



Matrix™ Diluent-2 (LISS)

Intended Use: _____

Matrix™ Diluent-2 (LISS) is Low ionic strength solution used for preparation of red cell suspension used for Matrix™ gel system.

Matrix™ Diluent-2 (LISS) kit components:

Matrix™ Diluent-2 (LISS)	Buffered low ionic strength solution of appropriate sodium chloride molarity.
Other Accessories	Package Insert

Performance Characteristics:

Standardization: Follows FDA guidelines

FEATURES	BENEFITS
Convenient	Ready to use Appropriate sodium chloride concentration to ensure that there is no hemolysis
Single Reagent Pack	With a shelf life of twelve months Is adequately buffered
Enhances Sensitivity	It enhances sensitivity by potentiating Ag-Ab reaction

Storage / Stability	Temperature	Duration
Unopened kit	2-8°C	12 Months
Opened bottle	2-8°C	12 Months

Available Pack Sizes	
100 ml	250 ml
500 ml	