



Matrix™ Direct Antiglobulin Test Card

Intended Use: _____

Matrix™ Direct Antiglobulin Test card is used to determine, if the red cells has been coated in vivo with immunoglobulin IgG, complement components C3d or both.

Matrix™ Direct Antiglobulin Test Card kit components:

Matrix™ Direct Antiglobulin Test Card	Matrix™ Direct Antiglobulin Test card contains six microtubes, prefilled with a gel in a suitable buffer containing Anti-Human IgG and Monoclonal Anti-C3d and neutral gel in appropriate microtubes
Configuration	IgG-C3d-Ctrl- IgG-C3d-Ctrl / Two profile per card
Other Accessories	Card fitment and Package Insert

Performance Characteristics:

Standardization: Follows FDA guidelines

FEATURES	BENEFITS
Stable End Reactions	Read results hours later Can be saved for QA or peer review Can be photocopied for documentation
Objective Interpretation of results	Easy Interpretation of weak reactions Relatively skill free
Differential detection	Detects and differentiates between IgG & C3d activity
Small sample size	More flexible More convenient
Indigenously Manufactured	Easy availability Long shelf life Another first from Tulip
Tulip's Technical Support	Dependability

Storage / Stability	Temperature	Duration
Unopened kit	4-25°C	12 Months
Opened card	---	Use immediately

Available Pack Sizes
24 Cards