



## Matrix™ Erygen-ABO Red Blood Cell panel for reverse grouping

**Intended Use:** \_\_\_\_\_

Matrix™ Erygen-ABO comprises of a 0.8% reagent red blood cells of A1, B and O suspended in an isotonic medium containing a preservative which is employed in the confirmation of the results of forward grouping.

**Matrix™ Erygen-ABO Kit components:**

<b>Matrix™ Erygen-ABO - A<sub>1</sub> cells</b>	0.8% reagent red blood cells of Group A <sub>1</sub>
<b>Matrix™ Erygen-ABO - B cells</b>	0.8% reagent red blood cells of Group B
<b>Matrix™ Erygen-ABO - O cells</b>	0.8% reagent red blood cells of Group O
<b>Other Accessories</b>	Package Insert

**Performance Characteristics:**

<b>Specificity:</b> ABO isoagglutinins	<b>Standardization:</b> Follows FDA guidelines
--	--

FEATURES	BENEFITS
<b>Ready to use</b>	Solution comprises of RBCs in an appropriate concentration in a preservative solution Components in precise quantities ensure maximum benefits during blood grouping
<b>Four component pack</b>	Separate reagents for different cell types

Storage / Stability	Temperature	Duration
<b>Unopened kit</b>	2-8°C	36 days
<b>Opened kit</b>	2-8°C	36 days

Available Pack Sizes
3 X 2 ml