



### Matrix™ ERYGEN-AS Antibody screening cell panels for Matrix™ Gel Cards

**Intended Use:** \_\_\_\_\_

Matrix™ ERYGEN-AS is 0.8% ready to use reagent red blood cells for Antibody Screening on Matrix™ Gel Cards.

**Matrix™ ERYGEN-AS kit components:**

<b>Matrix™ ERYGEN-AS</b>	3 cell panels each containing 0.8% red blood cells for Antibody Screening on Matrix™ Gel Cards.
<b>Other Accessories</b>	Antigen Chart, Package Insert

FEATURES	BENEFITS
<b>Convenient</b>	Available in a hassle-free, easy to use and ready prepared solution Preservatives are in precise quantities ensure maximum shelf life of antigens
<b>Single reagent pack</b>	Having a shelf life of 35 days and with added preservatives to maintain cell integrity

Storage / Stability	Temperature	Duration
<b>Unopened kit</b>	2-8°C	35 days
<b>Opened kit</b>	2-8°C	35 days

Available Pack Sizes
3 x 2 ml