



Matrix™ ERYGEN-PO Pooled reagent red cells for Antibody screening on Matrix™ Gel Cards

Intended Use: _____

Matrix™ ERYGEN-PO is 0.8% Pooled 'O' reagent red blood cells for Antibody Screening on Matrix™ Gel Cards.

Matrix™ ERYGEN-PO kit components:

Matrix™ ERYGEN-PO	0.8% Pooled 'O' reagent red blood cells
Other Accessories	Antigen Chart, Package Insert

FEATURES	BENEFITS
Convenient	Available in a hassle-free, easy to use and ready prepared solution Preservatives are in precise quantities ensure maximum shelf life of antigens
Single reagent pack	Having a shelf life of 35 days and with added preservatives to maintain cell integrity

Storage / Stability	Temperature	Duration
Unopened kit	2-8°C	35 days
Opened kit	2-8°C	35 days

Available Pack Sizes
1 x 5 ml