



**MatrixTM
ABO / Rho (D) Forward and Reverse Grouping Card with Autocontrol**

Intended Use: _____

MatrixTM ABO / Rho (D) Forward and Reverse Grouping Card with Autocontrol is used for determination of ABO and Rho (D) blood group systems and reverse grouping for confirmation of ABO grouping. Contains monoclonal Anti-A, Anti-B and Anti-D (VI-).

MatrixTM ABO / Rho (D) Forward and Reverse Grouping Card with Autocontrol kit components:

MatrixTM ABO / Rho (D) Forward and Reverse Grouping Card with Autocontrol	Matrix TM ABO / Rho (D) Forward and Reverse grouping card with Autocontrol contains monoclonal Anti-A, Anti-B and Anti-D (IgM) (VI-) with gel matrix from microtube 1-3. Microtube 4 serves as negative control. Microtube 5 and 6 contains neutral gel used for reverse grouping. The gel contains a suitable preservative
Configuration	A-B-D-Ctrl-A ₁ -B / One profile per card
Other Accessories	Card fitment and Package Insert

Performance Characteristics:

Standardization: Follows FDA guidelines

FEATURES	BENEFITS
Stable End Reactions	Read results hours later Can be saved for QA or peer review Can be photocopied for documentation
Objective Interpretation of results	Easy Interpretation of weak reactions Relatively skill free
Enhanced Sensitivity	Better detectability especially in antibody screen and cross match
Small sample size	More flexible More convenient
Indigenously Manufactured	Easy availability Long shelf life Another first from Tulip
Tulip's Technical Support	Dependability

Storage / Stability	Temperature	Duration
Unopened kit	4-25°C	12 Months
Opened card	---	Use immediately

Available Pack Sizes
24 Cards/24 Profiles