



Matrix™ Octoplus Complete Grouping Card

Intended Use: _____

Matrix™ Octoplus Complete Grouping Card is used for determination of ABO and Rho (D) blood group systems with autocontrol and reverse grouping for confirmation of ABO grouping.

Matrix™ Octoplus Complete Grouping Card kit components:

Matrix™ Octoplus Complete Grouping Card	Matrix™ Octoplus Complete Grouping Card contains Monoclonal Anti-A, Anti-B, Anti-AB and Anti-D (VI-) pre-filled with gel from microtube 1-4, Microtube 5 contains neutral gel which serves as a negative control and Microtube 6, 7 & 8 contains neutral gel for reverse grouping
Configuration	A-B-AB-D-Ctrl-A1-B-O / One profile per card
Other Accessories	Card fitment and Package Insert

Performance Characteristics:

Standardization: Follows FDA guidelines

FEATURES	BENEFITS
8 Microtubes	Forward grouping along with Anti-AB Reverse grouping including extra microtube for 'O' cells Complete blood grouping
Stable End Reactions	Read results hours later Can be saved for QA or peer review Can be photocopied for documentation
Objective Interpretation of results	Easy Interpretation of weak reactions Relatively skill free
Small sample size	More flexible More convenient
Indigenously Manufactured	Easy availability Long shelf life Another first from Tulip
Tulip's Technical Support	Dependability

Storage / Stability	Temperature	Duration
Unopened kit	4- 25°C	12 Months
Opened card	---	Use immediately

Available Pack Sizes
24 Cards