



MatrixTM Octoplus Forward and Reverse Grouping Card with Sub Grouping

Intended Use: _____

MatrixTM Octoplus Forward and Reverse Grouping Card with Sub Grouping is used for determination of ABO and Rho (D) blood group systems with subgrouping and reverse grouping for confirmation of ABO grouping.

MatrixTM Octoplus Forward and Reverse Grouping Card with Sub Grouping kit components:

MatrixTM Octoplus Forward and Reverse Grouping Card with Sub Grouping	Matrix TM Octoplus Forward and Reverse Grouping Card with Sub Grouping contains Monoclonal Anti-A, Anti-B, Anti-D (VI-) prefilled with gel from microtube 1-3, Microtube 4 contains neutral gel which serves as a negative control, Microtube 5 contains Anti-A1 lectin, Microtube 6 contains Anti-H lectin and Microtube 7 & 8 contains neutral gel for reverse grouping
Configuration	A-B-D-Ctrl-A1-H-A1-B / One profile per card
Other Accessories	Card fitment and Package Insert

Performance Characteristics:

Standardization: Follows FDA guidelines
--

FEATURES	BENEFITS
8 Microtubes	Forward grouping along with Anti-A1 and Anti-H Reverse grouping along with Autocontrol
Stable End Reactions	Read results hours later Can be saved for QA or peer review Can be photocopied for documentation
Objective Interpretation of results	Easy Interpretation of weak reactions Relatively skill free
Small sample size	More flexible More convenient
Indigenously Manufactured	Easy availability Long shelf life Another first from Tulip
Tulip's Technical Support	Dependability

Storage / Stability	Temperature	Duration
Unopened kit	4-25°C	12 Months
Opened card	---	Use immediately

Available Pack Sizes
24 Cards