



MatrixTM QAS Level-I

Intended Use: _____

0.8% suspension of red blood cells sensitized with IgG immunoglobulin for Quality Assurance testing of Anti-Human globulin reagent and related procedures in MatrixTM Gel System.

MatrixTM QAS Level-I kit components:

Reagent	0.8+/- 0.1% Suspension of Coombs Control Cells
Other Accessories	Package Insert

FEATURES	BENEFITS
Convenient	Available in a hassle-free, easy to use and ready prepared solution Components in precise quantities ensure maximum benefits during quality control procedures
Single reagent pack	Having a shelf life of 35 days and with added preservatives to maintain cell integrity

Storage / Stability	Temperature	Duration
Unopened kit	2-8°C	35 days
Opened kit	2-8°C	35 days

Available Pack Size
2 ml