



MatrixTM QAS Level-II

Intended Use: _____

Standardized, lyophilized monoclonal Anti-D IgG in specified concentration for Quality Assurance of reagents and techniques related to MatrixTM Gel System and classical tube technique.

MatrixTM QAS Level-II kit components:

Reagent	0.5 ml Standardized Anti-D IgG (Lyophilized) and 10 ml 6% BSA
Other Accessories	Package Insert

FEATURES	BENEFITS
Standardized Anti-D IgG	Ideal for Quality Assurance of reagent and pipetting techniques
Lyophilized reagent	Highly stable

Storage / Stability	Temperature	Duration
Lyophilized	2-8°C	12 months
Reconstituted	2-8°C	One week

Available Pack Size
0.5 ml Anti-D and 10 ml 6% BSA