



MatrixTM Rh Phenotype Card with Anti-D

Intended Use: _____

MatrixTM Rh Phenotype Card with Anti-D is used for determination of antigen of Rh system. The Rh phenotype determination is important for previously transfused patients, patients with irregular antibodies, during pregnancy. Contains monoclonal Anti-C, Anti-c, Anti-D, Anti-E, Anti-e.

MatrixTM Rh Phenotype Card with Anti-D kit components:

MatrixTM Rh Phenotype Card with Anti-D	Matrix TM Rh Phenotype Card with Anti-D contains six microtubes, pre-filled with gel in a suitable buffer with monoclonal Anti-C, Anti-c, Anti-D (IgM) (VI-), Anti-E, Anti-e and neutral gel. The gel contains a suitable preservative
Configuration	C-c-D-E-e-Ctrl / One profile per card
Other Accessories	Card fitment and Package Insert

Performance Characteristics:

Standardization: Follows FDA guidelines
--

FEATURES	BENEFITS
Stable End Reactions	Read results hours later Can be saved for QA or peer review Can be photocopied for documentation
Objective Interpretation of results	Easy Interpretation of weak reactions Relatively skill free
Enhanced Sensitivity	Better detectability especially in antibody screen and cross match
Small sample size	More flexible More convenient
Indigenously Manufactured	Easy availability Long shelf life Another first from Tulip
Tulip's Technical Support	Dependability

Storage / Stability	Temperature	Duration
Unopened kit	4- 25°C	12 Months
Opened card	---	Use immediately

Available Pack Sizes
24 Cards/24 Profiles