



Matrix™ Rh Phenotype Card with Anti-K

Intended Use: _____

Matrix™ Rh Phenotype Card with Anti-K is used for determination of Rh Phenotypes along with K antigen of Kell system. The Rh phenotype determination is important for previously transfused patients, patients with irregular antibodies, during pregnancy.

Matrix™ Rh Phenotype Card with Anti-K kit components:

Matrix™ Rh Phenotype Card with Anti-K	Matrix™ Rh Phenotype Card with Anti-K contains six microtubes, prefilled with a gel in a suitable buffer containing Monoclonal Anti-C , Anti-c, Anti-E , Anti-e, Anti -K and neutral gel in appropriate microtubes.
Configuration	C-c-E-e-K-Ctrl / One profile per card
Other Accessories	Card fitment and Package Insert

Performance Characteristics:

Standardization: Follows FDA guidelines

FEATURES	BENEFITS
Stable End Reactions	Read results hours later Can be saved for QA or peer review Can be photocopied for documentation
Objective Interpretation of results	Easy Interpretation of weak reactions Relatively skill free
Enhanced Sensitivity	High titre monoclonal antibodies ensures detection of weaker and heterozygous antigens
Small sample size	More flexible More convenient
Indigenously Manufactured	Easy availability Long shelf life Another first from Tulip
Tulip's Technical Support	Dependability

Storage / Stability	Temperature	Duration
Unopened kit	4-25°C	12 Months
Opened card	---	Use immediately

Available Pack Sizes
24 Cards