



## ReaCell I II III

### Antibody screening three cell panel for Matrix Gel Cards

**Intended Use:** \_\_\_\_\_

ReaCell I II III is 0.8% ready to use reagent red blood cells for Antibody Screening on Matrix Gel Cards.

**ReaCell I II III Kit components:**

ReaCell I II III	3 cell panel each containing 0.8% red blood cells for Antibody Screening on Matrix Gel Cards.
Other Accessories	Antigen Chart, Package Insert

FEATURES	BENEFITS
<b>Convenient</b>	Available in a hassle-free, easy to use and ready prepared solution Preservatives are in precise quantities ensure maximum shelf life of antigens
<b>Single reagent pack</b>	Having a shelf life of 53 days from the date of manufacturing and with added preservatives to maintain cell integrity

Storage / Stability	Temperature	Duration
Unopened kit	2-8°C	53 days
Opened kit	2-8°C	53 days

Available Pack Sizes
3 x 5 ml