



Tuliss®
Low ionic salt solution

Intended Use: _____

Tuliss®, a low ionic salt solution which is helpful in detection of weak antibodies during cross match techniques, antibody screening and antibody identification.

Tuliss® Kit components:

Tuliss®	Buffered, low ionic salt solution of appropriate sodium chloride molarity
Other Accessories	Package Insert

Performance Characteristics:

Standardization: Follows AABB and FDA guidelines

FEATURES	BENEFITS
Ready to use	No manual preparation required

Storage / Stability	Temperature	Duration
Unopened kit	2-8°C	24 Months
Opened kit	2-8°C	24 Months

Available Pack Sizes
5 ml